

Discover laboratory products at Rowe Scientific Pty Ltd.



A 100% Australian owned company supplying scientific laboratories throughout Australia and overseas since 1987.



*Haier Laboratory Refrigerators & Freezers*

**Haier Biomedical**

**IR0226 Pharmacy Refrigerator, Underbench, 118L**



HYC-118



**LED Display**

Microprocessor control, digital display and temperature adjustment with an increment of 0.1 C



**Multiple Alarm Functions**

Capable of alerting failures due to high and low temperature, sensor error, power failure, low battery, door ajar and remote alarm



**Ergonomic Design**

810mm cabinet height makes it easy to place under test bench, saving storage space.



**Lockable Design**

Safety lock to prevent unauthorized access



**Better Cabinet Temperature Uniformity, Auto Defrost**

The test point of performance in the cabinet increased from 2 points to 9 points.



**Energy Efficient Hydrocarbon (HC) Refrigeration System**

Optimised refrigeration system, utilising high-efficiency HC compressor significantly reduces energy consumption.

**Refrigerators, Spark Free**

**IR0225  
118L**



HLR-118SF



**IR0251  
310L**




HLR-310SF



**360° Explosion Proof Interior**

with ATEX II C-T6 EU Explosion Proof Certification.



**Setting Temperature Range +3°C to +16°C**

with precise control, keeps internal temperature variation within 3°C for great temperature uniformity.



**Clear and Complete Alarm System with One Key Search**

multiple alarms and data function which allows users to check and save historical data, which is useful and convenient for laboratory inspection processes.

See more at [www.rowe.com.au](http://www.rowe.com.au)

# Haier Biomedical

## -30°C Freezers, Spark Free

### IF0383 92L



DW-25L92SF



#### Explosion Proof Interior

SF series freezers feature standard interior design, door frame, and drawer. FL series freezers feature antistatic interior design, door frame, and drawer. All metal parts are grounded reliably and interior design is explosion proof. Certificates include China anti-explosion proof and EU explosion proof ATEX II C-T6



#### HC System to save energy

- Imported hydrocarbon compressor, and hydrocarbon refrigerants complete a true green, environmentally safe refrigeration system
- Cyclo-isopentane polyurethane foam insulation meets EU ROHS directive for a green and safe environment
- Optimized hydrocarbon refrigeration system design ensures a more energy efficient operation



#### Excellent temperature uniformity

Cold shelf evaporator design allows a faster cool down, and more uniform temperature distribution, typically  $\pm 4^{\circ}\text{C}$



#### Low noise

Unique structural design and system emits only 40 dba in sound level, eliminating unhealthy high pitch noise

### IF0421 278L



DW-30L278SF

**For more information visit our website**

**[www.rowe.com.au](http://www.rowe.com.au)**



**www.rowe.com.au**

Online, by phone or face to face, we give you the choice.

New South Wales & ACT 02 9603 1205 rowensw@rowe.com.au

Queensland 07- 3376 9411 roweqld@rowe.com.au

South Australia & NT 08- 8186 0523 rowsa@rowe.com.au

Victoria & Tasmania 03- 9701 7077 rowevic@rowe.com.au

Western Australia 08 - 9302 1911 rowewa@rowe.com.au

