



ROWE SCIENTIFIC

PTY LTD ABN 63 009 437 790

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987. We are proudly a 100% Australian owned company.

www.rowe.com.au



LIQUID DISPENSERS

Introducing the Dispensette® S by BRAND

The Dispensette® S range of liquid dispensers, with recirculating valve technology, are ideal for dispensing aggressive and non aggressive reagents safely, reliably and quickly.

FEATURES INCLUDE

- **Easy calibration technique**
Mechanical volume setting is easy to read; ensures accurate, repeatable volume control
- **Recirculating valve and tube**
Recirculation during priming avoids reagent waste
- **Large viewport**
Better control of liquid discharge and air bubbles
- **Hinged safety screw cap**
Doesn't get in the way when dispensing
- **Washer free sealing system**
Ease of cleaning and maintenance
- **Volume selection**
Volume setting is quick, secure and repeatable
- **360 degree rotating valve block**
Ensures bottle label is always facing the user
- **Telescoping filling tube**
Adjusts to different sized bottles.



FEATURES
RECIRCULATING
VALVE TECHNOLOGY

South Australia & Northern Territory
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

