



ROWE SCIENTIFIC

PTY LTD ABN 63 009 437 790

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987. We are proudly a 100% Australian owned company.

www.rowe.com.au



Quality Endorsed Company
ISO 9001:2015
UK 12512
SAS 0004

South Australia & Northern Territory
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au



LABORATORY HOTPLATES

Turning up the heat with laboratory hot plates

FEATURES

- Durable, thick cast aluminium top plate allowing for excellent heat transfer & uniform temperatures
- Cast top comes complete with drip edge for corrosion protection
- 316 stainless steel casing for corrosion resistance
- Durable rubber feet
- Separate control box fitted with simmerstat controller and heating light. Made from anodised aluminium, 2.5m rubber lead and plug for remote positioning
- Hotplate supplied with 2.5m, 10 Amp moulded mains lead and socket.

TECHNICAL SPECIFICATIONS

- Volts: 220-240V
- Hz: 50-60
- Power: 2.2Kw, single phase, 10 Amps
- Temperature: up to 350°C
- Heating Surface: 315mm x 585mm
- Overall height: 160mm



ROWE CODE: IH0062

Contact us now to find out more information about this product

www.rowe.com.au