



# Rowe Scientific PTY LTD

*For accuracy and professionalism*



**ADELAIDE**  
08 8186 0523  
rowesa@rowe.com.au

**BRISBANE**  
07 3376 9411  
roweqld@rowe.com.au

**MELBOURNE**  
03 8795 7771  
rowevic@rowe.com.au

**PERTH**  
08 9302 1911  
rowewa@rowe.com.au

**SYDNEY**  
02 9482 8433  
rowensw@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## QUANTUMCLEAN (Concentrate)

### Avoid Contamination Totally Clean Laboratory Glassware and Apparatus

For total cleaning (including radioactive decontamination) of all laboratory and industrial glassware, plastics, ceramics and stainless steel.



**IDEAL FOR USE IN ULTRASONIC CLEANERS  
PHOSPHATE FREE AND BIODEGRADABLE**

**USE:**

Supplied as a liquid concentrate for dilution with water (approx. 2-5% is adequate for most laboratory applications). Ideal for manual cleaning or for use in an ultrasonic cleaner.

**ORDERING INFORMATION:**

Rowe Code	Pack Size	Price
CQ1000	1L	\$ 17.58
CQ1050	5L	\$ 72.07
CQ1100	20L	\$191.63

**PLEASE ASK FOR  
A FREE SAMPLE!**



Proudly manufactured in Australia by Rowe Scientific Pty Ltd.