




www.rowe.com.au

Ultra Low Retention filter tips racked and sterile. ZEROTIP® JET BIOFIL® 2 STERILE R


• Economic order quantity is 2 boxes

VOLUME/SIZE CODE Pack size/Colour Code


10uL DT0279 10 x 96 tips
RACKED ZEROTIP
(31.6mm LONG) Red Rack




10uL XL DT0280 10 x 96 tips
RACKED ZEROTIP®
(44.5mm LONG) Red Rack




20uL XL DT0289 10 x 96 tips
RACKED ZEROTIP®
(50.46mm LONG) Yellow Rack




100uL DT0281 10 x 96 tips
RACKED ZEROTIP®
(50.5mm LONG) Yellow Rack




200uL DT0286 10 x 96 tips
RACKED ZEROTIP®
(59.2mm LONG) Yellow Rack




300uL DT0290 10 x 96 tips
RACKED ZEROTIP®
(59.2mm LONG) Yellow Rack



1000uL DT0282 10 x 96 tips
RACKED ZEROTIP®
(86.3mm LONG) Blue Rack



1000uL XL (1250uL) DT0283 10 x 96 tips
RACKED ZEROTIP®
(105.1mm LONG) Blue Rack



Colour coded boxes make it easy to identify tip volumes at a glance.

South Australia & Northern Territory
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

For accuracy and professionalism



Also available in Standard filter tips.

Standard Filter tip pipette sizes (Sterile)	Rowe Code	Price per box (960)
10uL RACKED (31.6mm LONG)	DT0272	10 x 96 tips Red Rack
10uL XL RACKED (44.5mm LONG)	DT0273	10 x 96 tips Red Rack
20uL RACKED (50.5mm LONG)	DT0274	10 x 96 tips Yellow Rack
100uL RACKED (50.5mm LONG)	DT0275	10 x 96 tips Yellow Rack
200uL RACKED (59.2mm LONG)	DT0276	10 x 96 tips Yellow Rack
300uL RACKED (59.2mm LONG)	DT0288	10 x 96 tips Yellow Rack
1000uL RACKED (86.3mm LONG)	DT0277	10 x 96 tips Blue Rack
1000uL XL (1250uL) RACKED (105.1mm LONG)	DT0278	10 x 96 tips Blue Rack



www.rowe.com.au



Ultra Low Retention filter tips racked and sterile

Discover the benefits of Biofil JET ZEROTIP® filter tips. Stock NOW available. Contact your local Rowe Scientific representative to request a sample.

- The surfaces of JET ZEROTIP® tips are produced through a special process, making the tip inner surface become super hydrophobic
- All tips are manufactured from super clear high quality Polypropylene
- Suitable for most brand micropipettors
- Colour-coded racks for easy tip volume identification
- Every rack or case is printed with lot No. for quality traceability
- DNase/RNase-free and non-pyrogenic
- Sterilized by irradiation SAL10-6 (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Using JET ZEROTIP® tips significantly reduces sample loss, and provides substantially higher reproducibility when working with critical media.

Also available:

Is a full range of single use plasticware for use in tissue culture and molecular biology including flasks and plates, 15 and 50mL centrifuge tubes and PCR tubes, serological pipettes plus cell strainers and inserts. Ask your LOCAL Rowe representative for details.

Discover more at www.rowe.com.au

South Australia & Northern Territory
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

For accuracy and professionalism

