



ROWE SCIENTIFIC
PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



Australian made surgical face masks.



ROWE CODE
HM0196
HM0042



These Australian made surgical face masks are manufactured with a 3 layer design along with an ultra-soft non-reactive inner layer to maximize facial comfort.

Manufactured to Australian standard: AS 4381:2002


| CODE | DESCRIPTION | COLOUR | UOM |
|--------|-------------------------------------------------|--------|--------|
| HM0196 | Mask, Surgical Face, Fluid Resistant With Ties | Green | 6 x 50 |
| HM0042 | Mask, Surgical Face, Fluid Resistant With Loops | Green | 6 x 50 |

Bacterial protection for both the wearer and the patient.

Med-con masks test at a PFE (particulate filtration efficiency) of 99.1% with a particle size of 0.1µm.

Breathability and comfort whilst wearing the mask is ensured with less than 1.5 units of delta pressure differential.

May also be used in an industrial application where similar particulate bacterial protection is required.



Support quality Australian made products.

Made to Australian standards, these masks are guaranteed to meet the highest level of quality to the end user.

The automated process of manufacturing reduces potential bio-burden contamination therefore ensuring consistency and guaranteed batch-to-batch regularity.

Material selection combined with Ultrasonic Sealing provides long time comfort and performance for the user.



These face masks may be used for source control, to help slow the spread of infection when a person speaks, coughs, or sneezes. Health care facilities should be aware that the use of face masks is different from personal protective equipment for health care personnel.

www.rowe.com.au

South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au