

24-Hour Emergency Number
(ChemTel) 888.255.3924 in U.S.

cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices



TERGAJET® Low-Foaming Powdered Detergent

- Concentrated to save you money
- Phosphate-free, biodegradable and readily disposable
- For FDA cleaning validations and food contact surface cleaning
- Low-foaming heavy duty detergent for high performance results
- Free rinsing to provide reliable results and no interfering residues
- Ideal for under-counter and floor-standing lab washing machines

Used to clean: Water testing labware, phosphate-sensitive lab equipment, hospital instruments, optical components, pharmaceutical and cosmetics manufacturing equipment, tissue culture ware, plate glass, metal castings, forgings and stampings, industrial parts, animal cages, electronic components, tanks and reactors.

Used to remove: Soils, radioisotopes, radioactivity, organics, grit, grime, slime, grease, tars, resins, fats, oils, blood, tissue, particulates, inorganic residues, chemicals and solvents.

Surfaces cleaned: Recommended for glass, stainless steel, PTFE, ceramic, cement and alkali resistant plastics and metals. TERGAJET is inhibited for aluminum. Corrosion testing is advisable.

Cleaning method: Labware washers, dishwashers, washer-sterilizers, washer-sanitizers, cart washers, spray cabinet washers, electrocleaning, ultrasonic cleaners, clean-in-place and manual cleaning.

Directions: Follow machine manufacturers directions and use the powder detergent receptacles. In general, use 1/2–1% in hot water wash cycle solutions (1 1/4–2 1/2 tbsp per gal, 3/4–1 1/4 oz per gal, or 5–10 g per L). Most machines need 2 1/2 oz added at a typical 2 gal wash cycle. For manual, soak, circulate or ultrasonic use, make a 1% solution in warm or hot water. For abrasive scrubbing, make a paste.

For difficult soils, raise water temperature and use more detergent. Wear protective gloves and eyewear when handling. RINSE THOROUGHLY—

preferably with running water. For critical cleaning, do final or all rinsing in distilled, deionized or purified water. Used on a wide range of glass, ceramic, plastic and metal surfaces. Corrosion testing is advisable. Material compatibility guidelines available on request.

| Convenient Sizes | Alconox Cat. # |
|---|----------------|
| Case 9 x 4 lb (9 x 1.8 kg) boxes | 2204 |
| 25 lb (11 kg) carton | 2225 |
| 50 lb (22 kg) carton | 2250 |
| 100 lb (45 kg) drum | 2201 |
| 300 lb (136 kg) drum | 2203 |
| 1 lb makes 26 gal (1 kg makes 200L) cleaning solution | |



TERGAJET passes inhibitory residue test for water analysis and is FDA certified. TERGAJET is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox Inc. detergents, use the Locate Dealers search tool in the Support section of alconox.com. To request FREE samples, use the Get Samples form at alconox.com. Or call Alconox Inc. at +914.948.4040, today.

TJTB.2

24-Hour Emergency Number
(ChemTel) 888.255.3924 in U.S.
cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices

| Physical Data | Typical Value |
|--|--|
| pH of 1% Solution | 11.5 |
| Surface Tension 1% Sol'n (dyne/cm) | 33 |
| Conductivity 1% Sol'n (mS/cm at 21 degrees C) | 8.9 |
| Form | Powder |
| Color | White |
| Solubility in Water | Completely soluble to 10% (w/w) |
| Hard Water Effectiveness | Highly effective |
| Biodegradability | Yes, meets EC648/2004 |
| Foam Tendency | Low foaming |
| Shelf Life | 24 months from end of month of manufacture |
| Fragrance Content | 0% |
| Phosphate Content (as phosphorus) | 0% |
| Flash Point (degrees F) | None |
| Organic Carbon (% calculated w/w) | 8.5% |

Chemical Description

Tergajet is an alkaline detergent containing sodium carbonate, silicate of soda, and surfactants. Tergajet contains nonionic surfactants.

Cleaning Validation and Residue Detection

To quantify trace levels of Tergajet for cleaning validation purposes, use conductivity, TOC, or HPLC. For details see the Cleaning Validation page in Support at alconox.com. For lab accreditation tests to show the detergent is rinsed away, see the Accreditation Residue Detection in the Support section of alconox.com.

Health and Safety Information

| | |
|--|--|
| OSHA Hazardous Ingredients | None |
| RCRA Hazard Class | Nonhazardous |
| EPA Priority Pollutants | None |
| DOT Hazard Class | Nonhazardous |
| Flammability | Nonflammable |
| Latex Content | None in product or packaging. |
| Irritation | Cat. 2 skin, Cat. 2A eye irritant |
| Oral Toxicity | LD ₅₀ > 5000 mg/kg oral rat. |
| Inhalation Toxicity | Solution: No, Powder dust: irritant |
| VOC Content | 0% |
| GMO Content | None |
| BSE/TSE Source | None |
| 2007/68/EC Annex IIIa Allergens | None |
| REACH Art. 57 (SVHC) | None |
| RoHS | No RoHS hazardous ingredients. |
| Carcinogenicity | NTP = No IARC = No OSHA = No |
| Ingredients on TSCA Inventory | Yes |
| EC 907/2006 | Contains no enzymes, optical brighteners, perfumes, allergenic fragrances, or preservative agents. |

Precautions

Avoid contact with skin, eyes and clothing. Do not take internally. In case of contact immediately flush skin with plenty of water. For eyes, flush with plenty of water for at least 15 minutes and get medical attention immediately. A Safety Data Sheet is available at alconox.com.

Contact Alconox Inc. for purchase specifications. Information presented is typical. Not to be taken as specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox Inc. detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

TJTB.2