




# MF SERIES LABORATORY REFRIGERATOR


**NULINE**  
REFRIGERATION

## FEATURES

- Top Mounted System
- Cabinet finish in Stainless Steel Interior/Exterior
- Automatic Defrost and Water Disposal
- 5 Adjustable Shelves per Door
- Available in Glass Door (TNG) or Solid Door (TN)
- Digital Display with Built-in Alarms
- Hi/Lo Readout



	2°C to 8°C
	700 L 1400 L
	740w x 830d x 2050h 1480w x 830d x 2050h



# MF SERIES LABORATORY REFRIGERATOR

**NULINE**  
REFRIGERATION

MODEL	MF70TN	MF70TNG	MF140TN	MF140TNG
Gross Capacity	700 L	700 L	1400 L	1400 L
External Dimensions in mm (w/d/h including castors or legs)	740 x 830 x 2050	740 x 830 x 2050	1480 x 830 x 2050	1480 x 830 x 2050
Internal Dimensions in mm (w/d/h)	600 x 690 x 1475	600 x 690 x 1475	1340 x 690 x 1475	1340 x 690 x 1475
Packed Dimensions in mm (w/d/h including pallet)	760 x 850 x 2100	760 x 850 x 2100	1500 x 850 x 2100	1500 x 850 x 2100
Gross Weight – kg (packed weight)	129	134	215	225
Net Weight – kg	115	120	191	201
Voltage	240V	240V	240V	240V
Temperature Range	2°C to 8°C	2°C to 8°C	2°C to 8°C	2°C to 8°C
Ambient Temperature Range	10°C to 32°C	10°C to 32°C	10°C to 32°C	10°C to 32°C
Refrigerant	R134a	R134a	R134a	R134a
Cooling System / Defrost	Fan Forced	Fan Forced	Fan Forced	Fan Forced
Exterior Cabinet Finish	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Inner Liner Material	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Door	Solid	Heated Glass	Solid	Heated Glass
Door Hinges	RH	RH	LH / RH	LH / RH
Door Handle	Full Length	Full Length	Full Length	Full Length
Type of Controller	Digital	Digital	Digital	Digital
Temperature Display	External	External	External	External
Data Logger	Optional	Optional	Optional	Optional
Door Ajar Alarms	Optional	Optional	Optional	Optional
High and Low Readout	Visual / Audible	Visual / Audible	Visual / Audible	Visual / Audible
Remote Alarm Contacts (BMS Connectivity for Alarms only) (BMS Connectivity for Temperatures)	Optional Optional	Optional Optional	Optional Optional	Optional Optional
Port Entry	Optional	Optional	Optional	Optional
Battery Back Up for Controller	Optional	Optional	Optional	Optional
Power Requirement	10amp GPO	10amp GPO	10amp GPO	10amp GPO
Castors	Yes	Yes	Yes	Yes
Adjustable Shelves	5	5	10	10
Usable Shelf Area in mm (w/d)	530 x 652	530 x 652	530 x 652	530 x 652
Light	LED	LED	LED	LED



**South Australia & NT**  
Ph: (08) 8186 0523  
rowesa@rowe.com.au

**Queensland**  
Ph: (07) 3376 9411  
roweqld@rowe.com.au

**Victoria & Tasmania**  
Ph: (03) 9701 7077  
rowevic@rowe.com.au

**New South Wales**  
Ph: (02) 9603 1205  
rowensw@rowe.com.au

**Western Australia**  
Ph: (08) 9302 1911  
rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)



Quality  
Endorsed  
Company  
ISO 9001:2015  
UK 10372  
SAS Global