

ENVIRONMENTAL & PERSONAL SAFETY SOLUTIONS

ROWE SCIENTIFIC PTY. LTD. ABN 63 009 437 790

EDITION
2



100% Australian Owned Since 1987

ENVIRONMENTAL & SITE SAFETY SOLUTIONS

TABLE OF CONTENTS

PERSONAL SAFETY FROM HEAD TO TOE	2-12
Gloves - Hot & Cold Protection	2
Nitrile Disposable Gloves	3
Gloves	4
Reusable Gloves	5
Disposable Gloves	6
Sleeve & Shoe, Face & Head Protection	7
Laboratory Coats	7
Aprons, Coats & Gloves	8
Hazmat Suits	9
Respiratory Protection	10-11
Safety Glasses	12
CONSUMABLES	13-19
Autoclave Bags	13
Indicator Strips & Air Monitoring	14-15
Sharps & Hazardous Waste	16
Light & Heavy Duty Cleaning Wipes	17
Kimberley-Clark Professional Wipes	18-19
SITE SAFETY	20-31
Eye Wash Stations & Safety Showers	20-21
Bottle Carriers, Bench Protectors, Mats	22
Safety Waste Caps & Funnels	23-25
Safety Cabinets	26-27
Glove Boxes & Desiccators	28
Hazard Warning Labels & Signs	29
Hazard Clean Up & Spill Stations	30-31
HOW TO ORDER	32

EDITION 2

© COPYRIGHT ROWE SCIENTIFIC PTY LTD. 2025.

Editor: Ilse Kaissis Graphics and setting: Ilse Kaissis

DISCLAIMER: Although every effort is made to ensure items are correct at time of print, it is the clients' responsibility to check the exact specification of the item at the time of ordering. All items and information in this catalogue are subject to change without notice. No part of this catalogue can be copied in any way without the consent of Rowe Scientific Pty Ltd.

PERSONAL SAFETY FROM HEAD TO TOE

GLOVES FOR HOT & COLD PROTECTION

GLOVES - HOT & COLD PROTECTION		SIZE	ROWE CODE
CRYOKIT 400 - CRYOGENIC GLOVE (PAIR) <ul style="list-style-type: none"> Waterproof 400mm glove length 		7	HC0358
		8	HC0359
		9	HC0360
		10	HC0361
		11	HC0362
CRYOKIT 550 - CRYOGENIC GLOVE (PAIR) <ul style="list-style-type: none"> Waterproof 550mm glove length 		8	HC0363
		9	HC0364
		10	HC0365
		11	HC0366
CRYOLITE-HP - CRYOGENIC GLOVE (PAIR) <ul style="list-style-type: none"> -200°C to 250°C temperature range Anti-static / waterproof 380mm glove length 		8	HC0340
		9	HC0341
		10	HC0342
		11	HC0343
CRYOLITE-HPS - CRYOGENIC GLOVE (PAIR) <ul style="list-style-type: none"> -200°C to 250°C temperature range Heavy-duty / waterproof 380mm glove length 		8	HC0555
		9	HC0556
		10	HC0557
		11	HC0558
CRYOLITE-WP - CRYOGENIC GLOVE (PAIR) <ul style="list-style-type: none"> -200°C to 250°C temperature range Knitted cuff / waterproof 260mm glove length 		7	HC0559
		8	HC0551
		9	HC0552
		10	HC0553
		11	HC0554
CRYOPLUS 2.1 - CRYOGENIC GLOVE (PAIR) <ul style="list-style-type: none"> Waterproof Textured finish for better grip 	LENGTH		
	380mm	7	HC0367
	380mm	8	HC0368
	380mm	9	HC0369
	380mm	10	HC0370
	380mm	11	HC0371
	450mm	7	HC0372
	450mm	8	HC0373
	450mm	9	HC0374
	550mm	9	HC0378
	550mm	10	HC0379
	550mm	11	HC0380



CRYOKIT 400



CRYOLITE-HP



CRYOLITE-HPS



CRYOLITE-WP



CRYOPLUS 2.1

PERSONAL SAFETY FROM HEAD TO TOE

NITRILE DISPOSABLE GLOVES

ANSELL GLOVES	SIZE	QTY	ROWE CODE
NITRILE DISPOSABLE EXAMINATION GLOVES			
MICROFLEX - BLACK <ul style="list-style-type: none"> • Powder Free • 245mm Cuff Length • Fully Textured 	Small	100	HG0080
	Medium		HG0081
	Large		HG0082
	XL		HG0083
MICROFLEX - VIOLET BLUE <ul style="list-style-type: none"> • Powder Free • 295mm Extended Cuff Length • Textured Fingers • Latex Free 	Small	50	HG0166
	Medium		HG0167
	Large		HG0168
	XL		HG0180
MICROTOUCH NITRATLEX - BLUE <ul style="list-style-type: none"> • Powder Free • 250mm Extended Cuff Length • Textured Surface • Latex Free 	Small	100	HG18884
	Medium		HG18886
	Large		HG1889
	XL		HG1888
TOUCH N TUFF - BLUE <ul style="list-style-type: none"> • Powder Free • 240mm Cuff Length • Textured Fingers • Latex Free 	Small	100	HG2400
	Medium		HG2405
	Large		HG2410
	XL		HG2415
TOUCH N TUFF - GREEN <ul style="list-style-type: none"> • Powder Free • 300mm Extended Cuff Length • Textured Fingers • Latex Free 	Small	100	HG2040
	Medium		HG2041
	Large		HG2042
	XL		HG2043
TOUCH N TUFF - GREEN <ul style="list-style-type: none"> • Powder Free • 240mm Cuff Length • Smooth Surface • Latex Free 	Small	100	HG1891
	Medium		HG1893
	XL		HG1897
TOUCH N TUFF - GREEN <ul style="list-style-type: none"> • Lightly Powdered • 240mm Cuff Length • Textured Fingers • Latex Free 	Medium	100	HG1892
	Large		HG1894
	XL		HG1896
VERSATOUCH - BLUE <ul style="list-style-type: none"> • Powder Free • 240mm Cuff Length • Textured Fingers • Latex Free 	Small	100	HG0001
	Medium		HG0002
	Large		HG0003
	XL		HG0004



MICROFLEX - BLACK



MICROFLEX - VIOLET BLUE



MICROTOUCH NITRATLEX - BLUE



TOUCH N TUFF - GREEN



ROWESCIENTIFIC
PTY LTD www.rowe.com.au

PERSONAL SAFETY FROM HEAD TO TOE

GLOVES

ANSELL GLOVES		SIZE	QTY	ROWE CODE
LATEX DISPOSABLE GLOVES - NATURAL				
TOUCH N TUFF <ul style="list-style-type: none">• Powder Free• 240mm Cuff Length• Fully Textured	Small	100	HG0052	
	Medium		HG0053	
	Large		HG0054	
	XL		HG0055	
TOUCH N TUFF <ul style="list-style-type: none">• Powder Free• 240mm Cuff Length• Smooth Surface	Small	100	HG0047	
	Medium		HG0048	
	Large		HG0049	
	XL		HG0050	
NEOPRENE DISPOSABLE EXAMINATION GLOVES - GREEN				
MICRO-TOUCH AFFINITY <ul style="list-style-type: none">• Powder Free• 240mm Cuff Length• Textured Fingers• Latex Free	Small	100	HG2010	
	Medium		HG2015	
	Large		HG2020	
	XL		HG2025	



TOUCH N TUFF



MICRO-TOUCH AFFINITY

ANSELL GLOVES		SIZE	QTY	ROWE CODE
NITRILE CUT-PROTECTIVE CHEMICAL-RESISTANT GLOVES - GREEN				
ALPHATEC 58-735 <ul style="list-style-type: none">• Double Wall (High Visibility Liner) Level C Cut Protection• Chemical and Liquid Resistant• Ansell Grip• 350mm Cuff Length	6	Pair	HG0580	
	7		HG0279	
	8		HG0280	
	9		HG0281	
	10		HG0282	
	11		HG0283	
NATURAL LATEX RUBBER-COATED HEAVY-DUTY GLOVES - BLUE				
ALPHATEC 62-401 <ul style="list-style-type: none">• 100% Cotton Liner• Heat Resistant• Gauntlet Cuff• Crinkle Grip• 320mm Cuff Length	7	6 Pairs	HG0311	
	8		HG0278	
	9		HG0326	
	10		HG0327	
	11		HG0328	



ALPHATEC 58-735



ALPHATEC 62-401



ROWESCIENTIFIC
PTY LTD www.rowe.com.au

PERSONAL SAFETY FROM HEAD TO TOE

REUSABLE GLOVES

ANSELL GLOVES		SIZE	QTY	ROWE CODE
BUTYL POLYMER GLOVES				
ALPHATEC - BLACK • 350mm Cuff Length (Beaded) • Rough Finish • Latex Free	Size 7	Pair	HG2045	
	Size 8		HG1803	
	Size 9		HG1801	
	Size 10		HG1804	
	Size 11		HG0503	
BUTYL POLYMER / VITON GLOVES				
ALPHATEC - BLACK • 300mm Cuff Length (Beaded) • Smooth Finish • Latex Free	Size 8	Pair	HG1956	
	Size 9		HG1957	
	Size 10		HG1958	
FIBRE-GLASS / STAINLESS STEEL - INFUSED LINER - CUT PROTECTION				
HYFLEX 74-500 - BLUE • Knitted • Knitwrist Cuff • 225 - 270mm Cuff Length • Food Processing	Size 7	12 Pairs	HG0309	
	Size 8		HG0310	
	Size 9		HG0692	
	Size 10		HG0693	
KEVLAR KNITTED YELLOW LINER AND BLUE RUBBER COATING - MODERATE CUT PROTECTION				
POWERFLEX 80-600 • Knitwrist Cuff • Rubber Coated Palm Finish • 250-280mm Cuff Length • Heat Resistant 200°C	Size 7	Pair	HG0212	
	Size 8		HG0213	
	Size 9		HG0214	
	Size 10		HG0215	
	Size 11		HG0216	
HPPE LAMINATED FILM HIGHLY CHEMICAL-RESISTANT GLOVES				
ALPHATEC 02-100 - WHITE • 5 Layer Laminated • Gauntlet Cuff • Smooth Finish • 380mm Cuff Length • Latex Free	Size 8	12 Pairs	HG0570	
	Size 9		HG0571	
	Size 10		HG0572	
PREMIUM SILVER LINED LATEX DISPOSABLE GLOVES				
ALPHATEC - PINK • 290mm Cuff Length (Pinked) • Pebble Embossed • Chlorinated	Size 6.5	12 Pairs	HG1400	
	Size 7		HG1440	
	Size 7.5		HG1500	
	Size 8		HG1550	
	Size 8.5		HG1600	
	Size 9		HG1650	
	Size 9.5		HG1700	
	Size 10		HG1702	



ALPHATEC - BLACK



HYFLEX 74-500 - BLUE



ACTIVARMR 80-600



ALPHATEC 02-100 - WHITE



ALPHATEC - PINK

PERSONAL SAFETY FROM HEAD TO TOE

DISPOSABLE GLOVES

MEDIFLEX GLOVES	SIZE	QTY	ROWE CODE
NITRILE DISPOSABLE EXAMINATION GLOVES - BLUE			
NISENSE NITRILE GLOVES • Powder Free • Latex Free • 245mm Cuff Length • Textured Fingertips • NSPF	XS	100	HG0452
	Small		HG0087
	Medium		HG0088
	Large		HG0089
	XL		HG0453
NISENSE 300 NITRILE GLOVES • Powder Free • Latex Free • 300mm Long Cuff Length • Textured Fingertips • NSPF	Small	10 x 100	HG0478
	Medium		HG0479
	Large		HG0480
	XL		HG0481
	XXL		HG0534
NITRASOFT NITRILE BLUE GLOVES • Powder Free • Latex Free • 245mm Long Cuff Length • Textured Fingertips • Non Sterile • NISFT	XS	10 x 200	HG0731
	Small		HG0732
	Medium		HG0733
	Large		HG0734
	XL		HG0735
	XXL		HG0736
NITRILE BLUE GLOVES • Extra Thick • Powder Free • Latex Free • 245mm Cuff Length • Textured Fingertips • PFNB	XS	100	HG0446
	Small		HG0447
	Medium		HG0448
	Large		HG0449
	XL		HG0450
	XXL		HG0451
NITRILE DISPOSABLE EXAMINATION GLOVES - BLACK			
BLACK WOLF GLOVES • Softer Formulation • Powder Free • Latex Free • 245mm Cuff Length • Textured Fingertips	Small	100	HG0454
	Medium		HG0455
	Large		HG0456
	XL		HG0457
	XXL		HG0458
NITRILE DISPOSABLE EXAMINATION GLOVES - PURPLE			
NICLEAN 4.2 GLOVES • Free From Thiurams, Thiazoles and Carbamates for Sensitive Skin • Powder Free • Latex Free • 245mm Cuff Length • Textured Fingertips	XS	200	HG0415
	Small		HG0416
	Medium		HG0417
	Large		HG0418
	XL		HG0419
LATEX DISPOSABLE EXAMINATION GLOVES - WHITE			
FLEXI GLOVES • Chlorinated Inner Surface for Easy Donning • Powder Free • 245mm Cuff Length • Fully Textured	XS	100	HG0464
	Small		HG0465
	Medium		HG0466
	Large		HG0467
	XL		HG0468



NISENSE NITRILE



NITRILE BLUE



BLACK WOLF



NICLEAN 4.2



FLEXI



PERSONAL SAFETY FROM HEAD TO TOE

SLEEVE & SHOE, FACE & HEAD PROTECTION

SLEEVES & SHOE COVERS	SIZE	QTY	ROWE CODE
Sleeve arm protection, elastic cuff, clear PVC, 0.1mm thick, light duty	450mm	100	HA1050
Sleeve arm protection, PE, blue	400mm	2000	HS0166
Disposable Shoe Cover, overshoe, non-slip, blue	Universal	500	HS2050
Disposable Shoe Cover (Surefoot), PP, non-slip, blue	Universal	500	HS2009



ROWE CODE
HA1050

EAR PLUGS & MUFFS	QTY	ROWE CODE
Prosafe disposable ear plugs, uncorded, class 4	200	RE1005
Superfit disposable ear plugs, large	200	HE0050
Ear Muff Peltor H9B, neckband/behind head	EA	HE1405



ROWE CODE
HE1405

FACE, EYE & HEAD PROTECTION (EA)	ROWE CODE
DFS2 - Face shield, clear, over-glasses, anti-fog and blue protective coating, elastic headband, ARTG listed	HF0118
MSA Black Eagle - Face shield, clear cover with chin guard	HF0167
Brow guard with ratchet headband	AP1050
Bolle Sphere complete face shield with head gear and visor	NM0003
Bolle Sphere face shield, clear, replacement visor only	NM0004
Hard Hat, V6 vented, full brim , white, 6 point ratchet harness with single handed adjustment	NH0022
Hard Hat, V6 vented, white, 6 point pushlock harness adjustment, with earmuff slots	NH0023
Hard Hat, V9 vented, white, 6 point pushlock harness adjustment, with earmuff slots	NH0021



ROWE CODE
AP1050



ROWE CODE
NM0003

Laboratory Coats

LABORATORY COATS	SIZE	ROWE CODE
LABORATORY COAT • Polycotton (White) • Normal Cuff • 2 x Waist Pockets and 1 x Chest Pocket	XS	HC0852
	SMALL	HC0833
	MEDIUM	HC0834
	LARGE	HC0835
	XL	HC0836
	XXL	HC0837
LABORATORY COAT • Polycotton (White) • Ribbed Cuff • 2 x Waist Pockets and 1 x Chest Pocket	XS	HC0853
	SMALL	HC0838
	MEDIUM	HC0839
	LARGE	HC0840
	XL	HC0841
	XXL	HC0842



PERSONAL SAFETY FROM HEAD TO TOE

APRONS, COATS & GLOVES

APRONS	SIZE	QTY	ROWE CODE
Medium-duty apron, clear PVC , 1250mm long	MEDIUM	12	PA1920
Reinforced apron, heavy-duty, black PVC , 1145mm long *BEST SELLER*	ONE SIZE	EA	PA1925
White PVC coated NYLON bib apron, limited acid resistance, 1219mm long	914 x 1219mm	EA	PA1040
Reinforced PVC heavy-duty apron, black, without side ties	ONE SIZE	12	PA1926
MICROCHEM 3000 apron, yellow 214 with sleeves	SMALL	EA	NA0007
	MEDIUM	EA	NA0011
	LARGE	EA	NA0012
	XL	EA	NA0013
Supported NEOPRENE , black apron	900 x 1150mm	12	PA1908
T-Cryo Light CRYOGENIC apron - Adjustable clips and buckles, easy clip closure (other sizes available)	700 x 900mm	EA	HC0346
T-Cryo Plus CRYOGENIC apron - Adjustable clips and buckles, easy clip closure (other sizes available)	700 x 1200mm	EA	HC0352

See our Cryogenic glove range on page 2

ALUMINISED REFRACTORY COATS	SIZE	QTY	ROWE CODE
Aluminised ARAMID / heavyweight parrawool furnace smock, double breasted with tape fastening, full length (ankle)	LARGE	EA	AC1093
Aluminised ARAMID furnace coat with wool lining, knee length	2XL	EA	AC0015
ALUMINISED ARAMID apron with sleeves, 1100mm long	SMALL	EA	AP1006
	MEDIUM	EA	AP1007
	LARGE	EA	AP1008
Coat ALUMINISED kevlar preox	LARGE	EA	AC1095

ALUMINISED REFRACTORY GLOVES	LENGTH	QTY	ROWE CODE
Aluminised Aramid back, Aramid woven palm, kevlar	400mm	EA	HG1114



ROWE CODE
AP1006



ROWE CODE
PA1925



ROWE CODE
AC1093



ROWE CODE
HG1114

Log into **www.rowe.com.au** to see your special price

PERSONAL SAFETY FROM HEAD TO TOE

HAZMAT SUITS

SPC 4400 / SPC 4800 / SPC 4900 SPLASH TIGHT SUITS

LIGHT AND COMFORTABLE TO WEAR: Draeger SPC 4X00 liquid-tight splash suits offer reliable protection against the finest particles and powders, many concentrated inorganic acids and alkaline solutions, and are also resistant to a wide variety of organic liquid chemicals. They are CE-certified.

DRAEGER PROTECTIVE HAZMAT SUITS	SIZE	ROWE CODE
SPC 4400 Front opening overalls with hood, elasticated cuffs & ankles. Yellow CPM material. Suitable for simple tasks and single use	LARGE	NO0008
	XL	NO0009
SPC 4400 Front opening overalls with hood, elasticated cuffs & ankles. Orange CLF material. Suitable for simple tasks and single use	LARGE	NO0010
	XL	NO0011
SPC 4800 Back opening overalls with hood, face cuff and permanently integrated gloves and socks. Orange CLF material, type 3/4/5/6 chemical protection* suit	LARGE	NO0012
	XL	NO0013
SPC 4800 Antistatic back opening overalls with hood, face cuff and permanently integrated gloves and socks. Orange CLF material, type 3/4/5/6 chemical protection* suit	LARGE	NO0014
	XL	NO0015
SPC 4900 Full-body protective suit with visor, integrated gloves and socks. Ideal for use with internal compressed air breathing apparatus. Orange CLF material, type 3/4/5/6 chemical protection* suit	LARGE	NO0016
	XL	NO0017



ROWE CODE
NO0010



ROWE CODE
NO0014



ROWE CODE
NO0017

Dräger

CHEMICAL PROTECTION*

TYPE 3 Protects against fluids and is used for chemical protection and oil spill clean-ups.

TYPE 4 Protects against spray and is used for chemical protection and oil spill clean-ups.

TYPE 5 Protects against particles and is used by asbestos removalists and in wood and plastic processing.

TYPE 6 Provides limited protection against spray and is used by asbestos removalists and in wood and plastic processing.

COVERALLS	SIZE	QTY	ROWE CODE
Safeguard GP style Coverall, type 5/6 protection, white	XL	EA	HC0742
SMS coveralls, elastic cuffs, hood and ankles, Type 5/6 protection, white	MEDIUM	50	HC0637
	LARGE		HC0635
	XL		HC0636
	XXXL		HC0667
Microgard 2000 coverall, disposable with anti-static treatment, white	LARGE	40	HC0436
	XXL		HC0438



ROWE CODE
HC0436

PERSONAL SAFETY FROM HEAD TO TOE

RESPIRATORY PROTECTION

DRAEGER RESPIRATORY PROTECTION		PACK SIZE	ROWE CODE
ODOUR AND PARTICLE FILTERING DISPOSABLE DUST MASKS			
X-PLORE 1720 - Valved, FFP2, low breathing resistance, light blue, single use	10		NM0031
X-PLORE 1920 - Valved, FFP2, low fogging, single use, Small	10		NR0023
X-PLORE 1920 - Valved, FFP2, low fogging, single use, Medium/Large	10		NR0022
TWIN FILTER HALF FACE MASK			
X-PLORE 3300 - Comfortable wearing for long lasting protection, advanced headstrap design allows for use under a helmet, requires 1 pair X-plore bayonet filters	S		NX0003
	M		NX0002
	L		NM0011
X-PLORE 3500 - Comfortable wearing for long lasting protection, advanced headstrap design allows for use under a helmet, hypoallergenic, requires 1 pair X-plore bayonet filters *PREMIUM*	S		NM0032
	M		NM0022
	L		NM0033
BAYONET FILTERS - FOR USE WITH X-PLORE 3300/3500/5500			
P3 R PARTICLE FILTER	11 Pairs		NF0008
PURE P2 R PARTICLE FILTER - Protects against particles in plastering, cleaning and varnishing	Pair		NR0024
PURE P3 R PARTICLE FILTER	10 Pairs		NF0007
A2B2 GAS FILTER - Organic and inorganic gases	9 Pairs		NF0006
A1B1E1 GAS FILTER - Respirator cartridge/filter	Pair		NF0022
A1B1E1K1 GAS FILTER	10 Pairs		NF0003
A2P3 R D COMBINATION FILTER	Pair		NR0025
A2B2P3 R D COMBINATION FILTER	6 Pairs		NF0005
A2B2E2K2Hg P3 R D COMBINATION FILTER - Respirator cartridge/filter	Pair		NF0021
A1B1E1K1Hg P3 R D COMBINATION FILTER - Respirator cartridge/filter	Pair		NR0011
FULL FACE MASK			
X-PLORE 5500 - Low fogging, 6-point harness, (mask only), requires bayonet filters (pair)	EA		NR0016
X-PLORE 6300 - Low fogging, 6-point harness, (mask only), requires RD40 filter	EA		NM0026
Carry Case for full face mask - Suits X-plore 5500/6300, hard case with shoulder strap	EA		NC0009
RD40 FILTERS - FOR USE WITH X-PLORE 6300 FULL FACE MASK			
1140 P3 R PARTICLE FILTER	EA		HP0140
1140 A2 P3 R D COMBINATION FILTER	EA		NR0026
1140 A2B2E2K1Hg P3 R D COMBINATION FILTER	EA		NR0027
1140 A2B2E2K2 HG P3 R D COMBINATION FILTER	EA		HR4043



ROWE CODE
NR0022/NR0023



ROWE CODE
NM0011/ NM0033



ROWE CODE
NR0024/NR0025/NR0011



ROWE CODE
NR0016



ROWE CODE
NM0026



ROWE CODE
HP0140

PERSONAL SAFETY FROM HEAD TO TOE

RESPIRATORY PROTECTION

DRAEGER RESPIRATORY PROTECTION		PACK SIZE	ROWE CODE
POWERED AIR PURIFYING RESPIRATOR (PAPR)			
PREMIUM KIT - Includes: X-plore 8500, blower, high capacity battery, belts and charger	Kit		NP0111
X-PLORE 8000 PAPR HELMET - With PC visor	EA		NH0020
X-PLORE 8000 PAPR FACE SHIELD	EA		NF0020
X-PLORE 8000 PAPR WELDING VISOR - Premium	EA		NW0001
FLEXIBLE BREATHING HOSE - Suitable for all headpieces	EA		NR0028
REUSABLE PARTICLE FILTER - Fits X-plore 8000 PAPR	EA		NR0029
Chemical filters in 6 different configurations - Call us!			
GAS DETECTION			
X-AM 2500 CSE KIT - Includes gas detector, pump and calibration gas	Kit		NG0202
PAC 6500 - O ₂ - Single gas detector for personal air monitoring (oxygen)	EA		NG0203
PAC 6500 - H ₂ S - Single gas detector for personal air monitoring (hydrogen sulfide)	EA		NG0204
PAC 6500 - CO - Single gas detector for personal air monitoring (carbon monoxide)	EA		NG0205
PAC 6500 - SO ₂ - Single gas detector for personal air monitoring (sulfur dioxide)	EA		NG0206
PAC 8000 - HCN - Single gas detector for personal air monitoring (hydrogen cyanide)	EA		IG0162
FILTERING ESCAPE DEVICE - RESPIRATOR			
PARAT 1260 - 10 minutes of protection using a multi-gas ABEK filter for organic and inorganic gases	EA		NE0014
PARAT 4720 - Industrial escape hood in soft case, 15 minutes of protection using a multi-gas ABEK P3 filter for dust, organic and inorganic gases	EA		NR0030
SELF CONTAINED BREATHING APPARATUS (SCBA)			
PSS 3000 - With 6.8L CF cylinder	1		NS0029



ROWE CODE
NH0020



ROWE CODE
NF0020



ROWE CODE
NW0001



ROWE CODE
NR0028



ROWE CODE
NR0029

Welding Protection

To comply with Australian standards when welding, you need to be wearing a Positive Air Purifying Respirator (PAPR). Your safety is delivered by Draeger through the combination of the X-plore 8500 blower with the X-plore 8000 Premium Welding Visor.

All that you need to work safely is combined in the Draeger Premium Welding Starter Kit, a bundle we have put together which comprises:

DRAEGER PREMIUM WELDING STARTER KIT ROWE CODE ND0041

X-PLORE 8500 BLOWER UNIT (IP)
X-PLORE 8500 HIGH CAPACITY BATTERY
X-PLORE 8000 WELDING BELT
X-PLORE 8000 HOSE COVER - Spark proof
FLEXIBLE BREATHING HOSE - Helmet/Visor/Hood
X-PLORE 8000 STANDARD CHARGER
X-PLORE 8000 WELDING VISOR - Premium
X-PLORE 8000 REUSABLE PARTICLE FILTER



ROWE CODE ND0041

PERSONAL SAFETY FROM HEAD TO TOE

SAFETY GLASSES

BOLLE SAFETY GLASSES (EA)	LENS	ROWE CODE
Atom - Safety goggles, blue TPR frame, indirect vents top / bottom, Platinum coating anti-scratch / anti-fog	Clear	NG0100
Axis 2 - Anti-scratch / anti-fog	Smoke	NG0012
Bandido	Clear	PS1172
BL30 - Anti-scratch / anti-fog	Clear	NG0182
Blade	Clear	PS1176
Contour	Clear	PS1177
Contour - Dark gun frame, indoor / outdoor	Dark	NG0026
Contour - Black frame, Platinum lite, anti-scratch / anti-fog	Smoke	NG0024
Hurricane - Black frame, anti-scratch / anti-fog	Smoke	NG0030
Hurricane - Black frame	Silver flash	NG0035
Hurricane - Black frame	Blue flash	NG0032
Hustler - Gloss black frame, anti-scratch / anti-fog	Clear	NG0037
Iri-S Diopter - +2.5 lenses with anti-scratch / anti-fog, black/grey frame	Clear	NG0046
Overlight II - Anti-scratch / anti-fog	Clear	NG0050
Override - Over-the-glasses	Clear	PS0091
Pilot 2 - Safety goggles, PP / TPR frame, indirect vented, Platinum coating anti-scratch / anti-fog	Clear	NG0127
Prism - Anti-scratch / anti-fog	Clear	NG0055
Prism Seal - Anti-scratch / anti-fog (sealed)	Clear	NG0063
Prism Seal - Anti-scratch / anti-fog (sealed)	Smoke	NG0064
Rush - Bronze glasses for twilight use, anti-scratch / anti-fog	Bronze	NG0073
Rush - Anti-scratch / anti-fog	Clear	NG0071
Rush+	Clear	PS0230
Rush+ - Platinum coating anti-scratch / anti-fog	Smoke	NG0078
Rush+ - Contrast lenses, gradient bronze technology	Clear	NG0081
Rush+ Glow - Phosphorescent temples, platinum coating anti-scratch / anti-fog	Clear	NG0086
Rush+ Seal - Assembled with gasket and strap, platinum coating anti-scratch / anti-fog	Clear	NG0083
Rush+ Small - Pink / white frame, platinum coating anti-scratch / anti-fog	Smoke	NG0090
Rush+ Small - Black / blue frame, platinum coating anti-scratch / anti-fog	Clear	NG0087
Silex - Anti-scratch / anti-fog	Clear	NG0171
Silium+ - Platinum, anti-scratch / anti-fog	Clear	NG0093
Tryon OTG - Over-the-glasses, platinum coating anti-scratch / anti-fog	Clear	NG0175
Ultim8 - Sealed glasses with strap, platinum coating anti-scratch / anti-fog	Clear	NG0190



**ROWE CODE
NG0032**



**ROWE CODE
NG0050**



**ROWE CODE
NG0127**



**ROWE CODE
NG0064**



**ROWE CODE
NG0073**



**ROWE CODE
NG0083**

CONSUMABLES

AUTOCLAVE BAGS

AUTOCLAVE BAGS - PAPER DIMENSIONS (L x W x H)	QTY	SACHEL NUMBER	ROWE CODE
GUARDIAN PAPER STEAM STERILISATION BAGS			
195 x 128 x 50mm	1000	1	PA1183
340 x 200 x 30mm	500	2	PA1175
340 x 250 x 60mm	500	4	PA1179
270 x 70 x 35mm	1000	7	PA1301
490 x 70 x 35mm	1000	9	PA1387
340 x 165 x 70mm	1000	13	PA1176
340 x 300 x 64mm	500	22	PA1173
440 x 202 x 98mm	250	30	PA1188
STERILOPE PAPER STEAM STERILISATION BAGS			
340 x 200 x 33mm	500	2	HA0003
340 x 250 x 60mm	500	4	HA0051
270 x 168 x 68mm	1000	5	PA1181
270 x 70 x 35mm	2000	7	PA0021
195 x 70 x 35mm	2000	8	PA9065
490x70x 35mm	1000	9	HA0047
155 x 90 x 50mm	2000	10	PA1186
340 x 165 x 70mm	1000	13	PA0023
440 x 202 x 100mm	500	30	PA1174
STERICLIN PAPER LAMINATE - LIGHT BLUE			
Pouch - Flat, transparent, 250 x 480mm	500		PP0118
Roll - Flat, see-through, 200mm x 200m	Roll		PS0637
Roll - Flat, transparent, 380mm x 200m	Roll		PS0352
Roll - Gusseted, transparent, 150mmW x 100m x 50mm gusset	Roll		PS0094
Roll - Gusseted, see-through, 250mmW x 100m x 65mm gusset	Roll		PS0588



ROWE CODE
PA1173



ROWE CODE
PA0021



ROWE CODE
PS0588

CONTAMINATED WASTE PLASTIC BAGS	QTY	ROWE CODE
Small 30.5 x 61cm, natural, printed, biohazard, plastic	500	PA1313
Small 30.5 x 60.5cm, natural, printed, biohazard	100	PA1320
Medium 57 x 75cm, 50 micron, natural, printed, biohazard	300	PA1314
Medium 57 x 75cm, 50 micron, YELLOW , printed, biohazard	300	PA0060
Medium 57 x 75cm, natural, printed, biohazard	100	PA1321
Large 72 x 91cm, 80 micron, natural, printed, biohazard	100	PA1315
Contaminated Waste Bag, 70cm x 180cm, 100 micron (RED PRINT)	50	PB1047
Waste Bag, biohazard, YELLOW , 54L, 32 x 90cmL, non-autoclavable	150	PB0326
Waste Bag, biohazard, YELLOW , 70cm x 112cm, autoclavable to 134°C	150	PB0272
Disposable, YELLOW , 380W x 940mmH, 120L	10	PW0095



ROWE CODE
PB0326

CONSUMABLES

INDICATOR STRIPS & AIR MONITORING

INDICATOR STRIPS & AIR MONITORING	QTY	ROWE CODE
TUBES		
Dräger Air Current Tube	10	HA1005
Ammonia 0.25/a, 0.25 - 3ppm	10	IC1790
Ammonia 5/a, 5 - 700ppm	10	IC1634
Ammonia 5/b, 5 - 100ppm	10	IC1702
Benzene 0.2 - 10ppm	10	NB0001
Carbon Dioxide 0.5 - 10% by volume	10	IC1611
Carbon Dioxide 1 - 20%	10	IC1613
Carbon Dioxide 0.5 - 6/0.1 - 1.2% by volume	10	IC1612
Carbon Dioxide 100 - 3000ppm (100/a)	10	IC1609
Carbon Dioxide 1 - 20%	10	IC1613
Carbon Dioxide 5 - 60% volume	10	IC1615
Carbon Disulfide 5/a, detector, 5 - 60ppm	10	IC1618
Carbon Monoxide 0.3 - 7% by volume	10	ND1610
Hydrochloric Acid 1/a, 1 - 10ppm	10	ID1679
Hydrogen Fluoride 1.5b, 1.5 - 15ppm & 6 - 60ppm	10	IC1782
Hydrogen Sulfide 0.2 - 5ppm	10	IC1640
Hydrogen Sulfide 0.2 - 7% by volume	10	IC1651
Hydrogen Sulfide 0.2/b, 0.2 - 6ppm	10	IC1645
Hydrogen Sulfide 2 - 60ppm	10	ND1647
Hydrogen Sulfide 0.5/a, 0.5 - 15ppm	10	IC1642
Hydrogen Sulfide 100/a, 100 - 2000ppm	10	IC1639
Hydrogen Sulfide 1/d, 1 - 200ppm / 1 - 20ppm	10	IC1674
Hydrogen Sulfide 2/a, 2 - 200ppm	10	IC1641
Mercaptan 20/a, 20 - 100ppm	10	IC1796
Mercury 0.1/b, 0.05-2mg/m3	10	ID1676
Natural Gas Test	10	IC1607
Water Vapour 0.1/a, 0.1 - 1.0mg/L	10	IC1603
GAS DETECTOR PUMP / ACCESSORY		
Accuro gas detection pump	EA	IC1690
Aspirator bulb for air-flow test tubes (eg smoke tube) (Dräger)	EA	ND0009

Dräger



ROWE CODE
IC1611



ROWE CODE
ND1610



ROWE CODE
IC1640



ROWE CODE
ID1676



ROWE CODE
IC1690

If you cannot see your required test, please call us

South Australia & NT
Ph: (08) 8186 0523

Queensland
Ph: (07) 3376 9411

Victoria & Tasmania
Ph: (03) 9701 7077

New South Wales
Ph: (02) 9603 1205

Western Australia
Ph: (08) 9302 1911

CONSUMABLES

INDICATOR STRIPS & AIR MONITORING

INDICATOR STRIPS & AIR MONITORING	QTY	ROWE CODE
CALIBRATION GAS		
CL2 10ppm in N2	58L	NG0005
NH3 75ppm in air	112L	NG0197
NH3 500ppm in air	112L	NG0159
NO2 10ppm in nitrogen	650L	NG0160
INDICATOR PAPER		
Litmus Blue	100	FI1400
Litmus Red	100	FI1450
pH 1-11 Tritest	5m	FI1735
pH 1-12 Duotest	5m	FI1700
INDICATOR STRIPS		
Ascorbic Acid 0-2000mg/L Quantofix	100	NI0005
Carbonate Hardness 0-20 degree Quantofix	100	NQ1060
Chloride 0-3000mg/L Quantofix	100	NQ1070
Formaldehyde 10-200mg/L Quantofix	100	NQ1200
Free Fatty Acids 0-3% Quantofix	50	NI0058
Glucose Semi-Quantitative 0-50-100-250-500-1000-2000mg/L Quantofix	100	NI0059
Iron 5-1000mg/L with reagents Quantofix	100	NQ1321
Medium Level Free Chlorine 10-200ppm	100	FI1829
Nitrate 10-500mg/L Nitrite 1-80mg/L Quantofix	100	NQ1530
Peracetic Acid 0-500mg/L Quantofix	100	NI0008
Peroxide 0-100mg/L Quantofix	100	NQ1600
Peroxide 0-25mg/L Quantofix	100	NQ1610
pH 0-14 in flat box pack (non bleeding) *BEST SELLER*	100	FI1100
pH 0-14 in snap cap tube *BEST SELLER*	100	FI1101
Phosphate 3-100mg/L with reagents Quantofix	100	NQ1620
Sulfite 0-1000mg/L Quantofix	100	NQ1720
Zinc 0-100mg/L with reagents Quantofix	100	NQ1960



ROWE CODE
NG0005



ROWE CODE
FI1400



ROWE CODE
FI1100



ROWE CODE
FI1101



ROWE CODE
NQ1720

***This is a brief selection from the available range.
Call us for the full range.***



ROWESCIENTIFIC
PTY LTD www.rowe.com.au

CONSUMABLES

SHARPS & HAZARDOUS WASTE

SHARPS CONTAINERS - FITTANK	VOLUME (L)	ROWE CODE
<ul style="list-style-type: none"> • Cytotoxic • Screw top • PURPLE • AS 4031:1992 compliant 	1.8 (Square)	PS0484
	8 (Oval)	PS0485
	19 (Square)	NS0030
<ul style="list-style-type: none"> • Screw top with petal insert and needle notch • YELLOW • AS 4031:1992 compliant 	1.8 (Square)	PS0486
<ul style="list-style-type: none"> • Screw top with petal insert • YELLOW • AS 4031:1992 	3 (Oval)	PC5001
	5 (Oval)	PC5002
	8 (Oval)	PC5003
	12 (Oval)	PC5004
	19 (Oval)	PS0616



ROWE CODE
PS0484



ROWE CODE
PS0486



ROWE CODE
PC5003

SHARPS CONTAINER	VOLUME (L)	ROWE CODE
<ul style="list-style-type: none"> • 340H x 265W x 195mmD • 50mm opening • Yellow • 2 piece cap (needle de-notching and screw cap) 	13	PC0554

BENCHTOP WASTE DISPOSAL	ROWE CODE
Benchtop Biohazard Disposal Can with Motion Sensor Lid - 1.5L, holds 25 x 25cm biohazard autoclavable disposal bag (LB0001 and NC0057)	NB0013
Biohazard Autoclavable Disposal Bag - Red with black biohazard symbol, HDPE, 0.72 mil thick, 22W x 28cmH, 100 pack	LB0001
Clavies Transparent Autoclavable Bags - 2 mil thick, 25 x 25cm, PP, 100 pack	NC0057



ROWE CODE
NB0013



ROWE CODE
LB0001



ROWE CODE
NC0057

CONSUMABLES

LIGHT & HEAVY DUTY CLEANING WIPES

Kleenex and Scott are trademarks of Kimberley-Clark Professional (K-CP).
All K-CP products are certified FSC (Forest Stewardship Council International).

HAND CLEANERS & WIPES	QTY	ROWE CODE
KIMTECH^{PRO} SCIENCE		
Kimwipes - Fine grade, low lint, 120 x 210mm, 280 wipes	1 box	FK1002
Kimwipes - Fine grade, low lint, 120 x 210mm, 280 wipes (bulk pack)	30 boxes	FK1008
OTHER CLEANING PRODUCTS		
Capri Interleaved Paper Hand Towel - Ultra slim fold, white, 24 x 23cm, 150 sheets	16 packs	HT0374
Capri Paper Hand Towel Roll Premium - 100m x 18cm, white	16 rolls	HT2042
SCOTT Optimum Interleaved Hand Towel - 30.5cm x 21cm, white, 150 sheets	16 packs	HT2031
SCOTT Small Paper Towel Roll - 100m x 18cm, white	16 rolls	HT0310
SCOTT Long Paper Towel Roll - 140m x 18cm, white	8 rolls	HT0039
Kleenex Optimum Interleaved Hand Towel - 30.5 x 24cm, white, 120 sheets	24 packs	HT2022
GLR Alcohol Wipes - 70% Isopropyl Alcohol / 30% DI Water, 42 x 14cm (75 wipes) *BEST SELLER*	EA	HW0028
TISSUES		
Tissues Tork - 2 ply, 195 x 200mm, 210 tissues	48 boxes	FT1023
Tissues Kleenex - 2 ply, 195 x 205mm, 100 tissues	48 boxes	FT1025
Tissues Kleenex - 2 ply, 195 x 205mm, 200 tissues	24 boxes	FT1020
MERCURY DECONTAMINATION		
Hand cleaner cream for MERCURY decontaminant	340g	NH0024
Mercon wipes for MERCURY decontamination *BEST SELLER*	70	CM0088
MERCURY DECONTAMINANT POWDER HgX	9kg	CM3772
MERCURY SPILL KIT - 3 x use, absorbent capacity: components for 3 x spills	10L pail	CM3921



ROWE CODE
FE2000



ROWE CODE
FK1016



ROWE CODE
HT0374



ROWE CODE
HW0028













For Hazard Cleanup - Pads and Spill Kits - See Page 30-31

CONSUMABLES

KIMBERLEY-CLARK PROFESSIONAL WIPES

WypAll Product Selector Guide

 Primary Task  Secondary Task	Product Name/ Descriptor	 Wipe Water & Liquids	 Wipe Oil/Grease	 Durability	 Quick Absorbing	 Chemical Solvents	 Low Lint	 Dust & Dirt	 Special Benefits	ROWE CODE
GENERALCLEAN™ Light to medium-duty general cleaning solutions	LIGHT CLEANING TOWELS L10	●●●●○						✓	Light-duty, single-use towel	FK1016
	MEDIUM CLEANING TOWELS L20	●●●●○						✓	Single-use towel for medium-sized jobs	FK1012
	HEAVY CLEANING TOWELS L30	●●●●○	●●●●○		✓			✓	Soft, absorbent towel for quick cleanups	FK1006
	CLEANING CLOTHS X50	●●●●○	●●●●○	✓	✓			✓	Absorbent cleaning cloth for light industrial use	MANY see next page
	MULTI-TASK CLEANING CLOTHS X60	●●●●○	●●●●○	✓	✓	✓	✓	✓	Versatile cloth; absorbent and durable; assorted rag replacement	HW0065 HW0073
	MICROFIBER CLEANING CLOTHS	●●●●+		✓	✓			✓	Color-coded; designed to deliver superior cleaning performance	HW0067
CRITICALCLEAN™ Help reduce contamination risk with cleaning, disinfecting and sanitizing solutions	REGULAR-DUTY COLOR-CODED WIPES	●●●●○	●●●●○	✓					Color-coded for cross-contamination control. Suitable for light wiping tasks	HW0070
	HEAVY-DUTY COLOR-CODED WIPES	●●●●○	●●●●○	✓					Color-coded for cross-contamination control. Suitable for heavy wiping tasks	HW0074 HW0069
	WETTASK™ CUSTOMIZABLE WET WIPING SYSTEM FOR BLEACH, DISINFECTANTS AND SANITIZERS	●●●●○	●●●●○	✓		✓	✓		Meltblown Basesheet with broad chemical compatibility	HW0071
POWERCLEAN™ Maximize performance while minimizing time, effort and product waste	MEDIUM DUTY CLOTHS X70	●●●●○	●●●●○	✓	✓	✓	✓	✓	Replacement for rags, great for low-lint wiping	HW0041 HW0075
	HEAVY DUTY CLOTHS X80	●●●●+	●●●●+	✓	✓	✓	✓	✓	Heavy-duty, extra-durable cleaning cloth	HW0077 HW0078
	EPIC™ WIPERS	●●●●○	●●●●+	✓		✓	✓		Absorb large amounts of oil, grease and ink	FE2000
	WETTASK™ CUSTOMIZABLE WET WIPING SYSTEM FOR CHEMICALS AND SOLVENTS	●●●●○	●●●●○	✓	✓	✓	✓		Highly absorbent Hydroknit basesheet with chemical compatibility	HW0071

CONSUMABLES

KIMBERLEY-CLARK PROFESSIONAL WIPES

WYPALL CLEANING WIPES	QTY	ROWE CODE
GENERAL CLEAN		
L10 - Heavy duty for light cleaning, centrefeed wipes, white, 300m x 20cm	4 rolls	FK1016
L20 - Heavy duty centrefeed wipes, perforated, white, 165m x 20cm	4 rolls	FK1012
L30 - Heavy duty, 42 x 39cm, blue, 200 wipes	1 box	FK1006
X50 - Rag on a roll, reinforced, white, 70m x 49cm	3 rolls	FR1120
X50 - Rag on a roll, reinforced, white, 70m x 24.5cm	4 rolls	FR1125
X50 - Rag on a roll, reinforced, blue, 70m x 49cm	3 rolls	FR1100
X50 - Rag on a roll, reinforced, blue, 70m x 24.5cm	4 rolls	FR1200
X50 - Reinforced, single sheet wipes, blue, 32.5 x 49cm, 250 wipes	1 box	HW2610
X50 - Reinforced, single sheet wipes, white, 24 x 42cm, 160 wipes	4 boxes	HW2630
X50 - Reinforced, single sheet wipes, white, 24 x 42cm, 75 wipes	8 boxes	HW0049
X60 - Multi-purpose wipes, jumbo roll, blue, 31.5 x 31cm, 1100 wipes	1 roll	HW0065
X60 - Multi-purpose wipes, small roll, blue, 24.5 x 42cm, 80 wipes	6 rolls	HW0073
MICROFIBRE CLOTHS - Wash and re-use, green, 40 x 40cm, 6 per pack	4 packs	HW0067
CRITICAL CLEAN		
REGULAR DUTY WIPES - Green, 60 x 30cm, 20 wipes	12 packs	HW0070
HEAVY DUTY WIPES - Yellow, 34 x 45cm, 106 wipes	6 packs	HW0074
HEAVY DUTY WIPES - Red, 34 x 45cm, 106 wipes	6 packs	HW0069
WETTASK - Dry wipes for your choice of wetting agent (ask Rowe Scientific to formulate), white, 30.4 x 15.2cm, 140 sheets per roll	6 rolls	HW0071
POWER CLEAN		
X70 - Medium duty cloth wipes, white, 31.5 x 42.5cm, 160 wipes	1 pack	HW0041
X70 - Medium duty cloth wipes, white, 31.5 x 42.5cm, 220 centrefeed roll wipes	4 rolls	HW0075
X80 - Heavy duty cloth wipes, white, 31 x 34.5cm, 540 perforated wipes per roll	1 roll	HW0077
X80 - Heavy duty cloth wipes, white, 31.5 x 42.5cm, 150 wipes	Case	HW0078
EPIC WIPERS - Brag box wipers, blue, 42 x 34.5cm, 250 wipes	1 pack	FE2000
WETTASK - Dry wipes for your choice of wetting agent (ask Rowe Scientific to formulate), white, 30.4 x 15.2cm, 140 sheets per roll	6 rolls	HW0071



**ROWE CODE
FK1016**



**ROWE CODE
HW0049**



**ROWE CODE
HW0067**



**ROWE CODE
HW0069**



**ROWE CODE
HW0071**

WYPALL
BRAND



ROWE SCIENTIFIC
PTY LTD www.rowe.com.au

SITE SAFETY

EYE WASH STATIONS & SAFETY SHOWERS

EYE WASH STATIONS

ROWE CODE

Eye Wash Station Rack - Plastic (without eye wash bottle, requires **PE1050**)

PE1100

Eye Wash Bottle - PE, 500mL ***BEST SELLER***

PE1050

Eyewash Saline Solution - Bottle with eye cap

NE0008

Eye wash single hand held drench hose with 1.5m hose

HE0052

Safety Shower and Eye Wash - No foot treadle

HS1005

Pedestal mounted eye and face wash nozzle with bowl - No foot treadle

HP0227



**ROWE CODE
PE1100 & PE1050**



**ROWE CODE
HS1005**

**ROWE CODE
HP0227**

PRATT
SAFETY SYSTEMS



Portable Gravity Fed Eye Wash Unit
ROWE CODE IE0226



Portable Gravity Fed Eye Wash Unit
ROWE CODE IE0347



**Pratt Water
Preservative Solution**
250mL (4 pack)
Suitable only for gravity fed units.
ROWE CODE CE0308

EMERGENCY TREATMENT

ROWE CODE

Hexafluorine Solution - LPMF 500mL portable eyewash bottle (commercial formulation for HF acid splashes - To be used under medical supervision)

NH0003

Eye wash Afterwash - For use after hexafluorine wash (to be used under medical supervision)

CA0719

Diphoterine Solution - 5L portable shower with green storage case

CD0031

Diphoterine Solution DAP LPMD - 500mL eye wash bottle

CD0622

Diphoterine Solution Micro DAP - 100mL spray

CM1009

Diphoterine Solution Mini DAP - 200mL aerosol for chemical burns

CM1003

Diphoterine Solution LMPE.DAS Mural Wall Mount Wash Station - Includes 2x 500mL eyewash bottles, 1x 200mL mini spray, 1x 200mL afterwash II

CD0029

Calcium Gluconate Gel - 50g tube (Must have for users of hydrofluoric acid) ***BEST SELLER***

CC2350



**ROWE CODE
CD0029**

SITE SAFETY

EYE WASH STATIONS & SAFETY SHOWERS

SAFETY EQUIPMENT	ROWE CODE
Test Kit for Safety Showers (with bin)	HT0341
Safety Shower test sock with handle and receptacle Test Kit	HS1018
Emergency Shower and Eye Wash Sign - 450 x 300mm, metal (other sizes available)	HS0819



ROWE CODE
HS0819



ROWE CODE
HT0341

ROWE CODE
HS1018

Safety Shower & Eye Wash Station

ROWE CODE **HS1007**

FEATURES INCLUDE:

- **Structural components made of 304 stainless steel**
For strength and durability
- **316 stainless stay-open valves**
For dependable operation and performance
- **Shower head made of ABS plastic or optional 316 stainless steel**
With a unique conical water pattern to gently deliver large volumes of water ensuring complete decontamination
- **Shower flow control regulator**
Limiting shower output to 90LPM
- **Interchangeable nozzles**
Converting from an eye wash to an eye-face wash



ROWE CODE
HS1007



Deluge Plastic Shower Head
As standard on HS1007



Deluge Stainless Steel
Shower Head
As standard on HC0758



Aerostream Triple
Nozzle (Aerated)
HP0091



Aerostream Single
Nozzle (Aerated)
HP0272

SITE SAFETY

BOTTLE CARRIERS, BENCH PROTECTORS, MATS

BOTTLE CARRIERS (EA)

ROWE CODE

Bottle Carrier - 1 x 4L bottle, 180 x 180mm ID (white plastic coated wire)	HB1060
5L Flask Carrier - 23cm square base, 45cm overall height	HC1280
Bottle Carrier - 2 x 2.5L winchester bottles, 140 x 140mm, (white plastic coated wire)	HB1050
Bottle Carrier - 2 x 4L bottles, 165 x 165mm ID each side (white plastic coated wire)	HB1070
Bottle Carrier - 4 place, 310 x 300 x 270mm, holds up to 120mm DIA bottles	PB0086
Bottle Carrier - 6 x 500mL bottles, 75 x 75mm ID each section (white plastic coated wire)	HB1044
Bottle Carrier - 6 x 1L bottles, 100 x 100mm ID each section (white plastic coated wire)	HB1045
Bottle Carrier - 6 x 1L bottles, 105 x 105mm ID each section (white plastic coated wire)	HC5091
Bottle Carrier - 6 place, 320 x 300 x 200mm, holds up to 95mm DIA bottles	PB0085
Bottle Carrier - 12 place, 80 x 80mm ID each section (white plastic coated wire)	HB1035




ROWE CODE
PB0085 & PB0086

BENCHTOP PAPER PROTECTOR SHEETS

QTY

ROWE CODE

Benchtop Protection Paper Sheets (Labtop) - 48 x 60cm PE coated (MN 210 PE) (similar to Benchkote)	50	FL1005
Benchtop Protection Paper - 48 x 60cm, PE coated (MN 210 PE)	100	FL1000
Benchtop Protection Paper Roll (Labtop) - 60cm x 50m, PE coated (MN 210 PE)	1 Roll	FL1025
Benchkote Benchtop Protector Sheets - 46 x 57cm	50	FB1000
Benchkote Benchtop Protector Roll - 46cm x 50m	1 Roll	FB1010
Benchtop Protection Paper Mat - 4 ply, plastic backed, 41.5cm x 91m	2 Rolls	HM1290



SAFETY MATS (EA)

ROWE CODE

Safety Matting - Ergonomic wet/dry, 900 x 900mm, with drainage holes (premium)	HM0022
Bevelled Edges Female - 900 x 75mm, yellow (to suit HM0022)	HM0023
Bevelled Edges Male - 900 x 75mm, yellow (to suit HM0022)	HM0024



SITE SAFETY

SAFETY WASTE CAPS & FUNNELS

REDUCE YOUR EXPOSURE TO SOLVENTS

SCAT SAFETY CAPS

If you use solvents in your workplace, SCAT provides the tools to assist you in controlling your exposure, whilst maintaining the integrity of your analysis. The range extends from HPLC solvent delivery and resultant waste collection, to safety funnels - all with air valves (carbon filled) to reduce fugitive emissions.

The Safety Standard for HPLC

Safety Caps stop solvent vapours from escaping and keep your eluents clean. They protect your health, the environment and the quality of your analysis. The mixing ratio of your solvents remains stable, as there is no evaporation.

GL45 is the most common thread, however other sizes are also available.



SAFETY CAP III FOR HPLC SOLVENT SUPPLY • 3 X CAPILLARY CONNECTIONS • GL 45 THREAD • PTFE				ROWE CODE
	AIR VALVE*	SIZE FITTING* for capillary connection	BLIND PLUG*	
SINGLE SAFETY CAP	1	3 x each colour (9)	2	NS0036
4 CAP STARTER KIT	4	4 x 3 x each colour (36)	8	NS0037

*(UNF1/4" 28G thread)



AIR VALVES FOR SAFETY CAPS			ROWE CODE
<ul style="list-style-type: none"> • UNF 1/4" 28G thread • Service life indicator • PTFE filter 	1 per pack		NA0053
	10 per pack		NA0054
	100 per pack		NA0055



ACCESSORIES INCLUDED!

Fittings and blind plugs in every size are included in the scope of delivery.
ID: 1.6mm, 2.3mm, 3.2mm

FOR TRACE LEVEL ANALYSIS OF PFAS

A different supply system manufactured from polypropylene is required.

SAFETY CAP III PFAS FOR HPLC SOLVENT SUPPLY • 3 X CAPILLARY CONNECTIONS • GL 45 THREAD • GREEN PP				ROWE CODE
	AIR VALVE*	SIZE FITTING* for capillary connection	BLIND PLUG*	
SINGLE SAFETY CAP	1	3 x each colour (9)	2	NS0038
4 CAP STARTER KIT	4	4 x 3 x each colour (36)	8	NS0039

*(UNF1/4" 28G thread)

PFAS AIR VALVES FOR SAFETY CAPS			ROWE CODE
<ul style="list-style-type: none"> • PFAS Analytic • Service life indicator • Green PP • UNF 1/4" 28G 	1 per pack		NA0056
	10 per pack		NA0057



THE No. 1 System
for Solvent Safety

SITE SAFETY

SAFETY WASTE CAPS & FUNNELS

To safely trap solvent residues, and protect you from harmful vapours, choose the MARCO safety funnel. The hinged lid and integrated safety seal stops harmful vapours from escaping. The removable sieve captures large particles and stirring bugs. Made of HDPE, the Marco Safety Funnel will give you years of use.

MARCO SAFETY FUNNEL	BOTTOM THREAD	ROWE CODE
• 140mm funnel diameter • Hinged lid & removable sieve (HDPE)	GL38	PS0631
	GL45	PF0174
	S58	NS0031
Funnel to fit Safety Waste Cap "LISA"	M 30x35	NS0043



ROWE CODE
PF0174



Safe Closure

- Comfortable handling: easy, quick opening and closing with "click"-closure

Hinged Lid with Sealing Ring

- For secure closure, when the funnel is not in use

Stable Hinge

- A double-screw hinge, constructed of rust-free stainless steel: for long life and stability

Dirt Sieve

- For collecting stirring magnets or large dirt particles. Laboratory bottles can also be placed here to dry off. The hole size has been optimised in order to improve the run-off of chemicals. **Both sieve and splash protection can be removed for cleaning.**

140 mm

Fine Mesh Sieve

- Supplied with plastic coarse mesh (Fine mesh stainless steel optional)

Various Screw-Threads

- Suitable for different sizes of collecting container

SITE SAFETY

SAFETY WASTE CAPS & FUNNELS

SCAT SAFETY WASTE CAPS

*Having secured the integrity of your analysis, what do you do with the waste?
These waste mixtures are highly hazardous to your health.*

CALL US TO SECURE YOUR SAFETY.

If you are collecting into a bottle, or into a jerry can, we have a waste cap to suit your needs.

SAFETY WASTE CAP FOR HPLC SOLVENTS	ROWE CODE
GL45 connection, 2 x capillary, 1 x tubing (scat)*	NS0040
S60/61 connection, 2 x capillaries, 1 x tubing, mechanical level indicator*	NS0041
S60/61 connection, 2 x capillaries, 1 x tubing, electronic level indicator*	NS0042



ROWE CODE
NS0040

If you need to collect waste from additional sources, choose the Safety Waste Cap "Lisa", which allows for the addition of a "Marco" Safety Funnel.

SAFETY WASTE CAP "LISA" FOR HPLC SOLVENTS	ROWE CODE
GL45 connection, 4 x capillary fittings (each 3 sizes), 3 x tubing, all blind plugs*	NS0045
S58x4 connection, 4 x capillary fittings (each 3 sizes), 2 x tubing, all blind plugs *	NS0052
SAFETY FUNNEL "MARCO" FOR LISA SAFETY CAP	ROWE CODE
M 30x35 thread (bottom) with hinged lid (140mm ID) and removable sieve, HDPE	NS0043



ROWE CODE
NS0052

ROWE CODE
NS0043

THREAD ADAPTORS FOR DIFFERENT CONTAINERS	ROWE CODE
*GL45 (male) to GL38 (female), polypropylene. To connect to some 5L "jerry cans" - check thread	PA0079
*S60/61 (male) to S58x4 (female), polypropylene. To connect to Australian 10 and 20L "plastic drums"	PA0082



SCAT EXHAUST FILTERS

The efficacy of all of these systems is based on the safe filtering/binding of vapour/gases vented back into the laboratory (your breathing space). SCAT have extensively developed their exhaust filter to minimise your exposure.

To satisfy the various waste compositions, temperature and flow rates, SCAT have developed their exhaust filter in 3 sizes.

EXHAUST FILTER FOR SAFE VENTING OF HPLC WASTE GASES	ROWE CODE
With change indicator small - 3 month life (each)	NE0017
With change indicator medium - 6 month life (each)	NE0018
With change indicator medium - 6 month life (2 pack)	NE0019
With change indicator large - 12 month life (each)	NE0020
With change indicator large - 12 month life (2 pack)	NE0021



SITE SAFETY SAFETY CABINETS

CORROSIVE SUBSTANCE

Indoor Storage Cabinets - **Metal**



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
30	770 x 515 x 465	1	1	HC3220
60	1065 x 515 x 465	1	2	HC3250
100	770 x 935 x 620	2	1	HC3260
160	1295 x 1115 x 500	2	2	HC3263
250	1825 x 1115 x 500	2	3	HC0690
350	1825 x 1440 x 585	2	3	HC0691
425	1825 x 1115 x 865	2	3	HC0692



ROWE CODE
HC0691

CORROSIVE SUBSTANCE

Indoor Storage Cabinets - **Polyethylene**



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
20	600 x 460 x 466	1	1	HC3225
40	770 x 600 x 466	1	1	HC3205
80	1150 x 600 x 466	1	2	HC0518
100	770 x 880 x 606	2	1	HC3207
110	770 x 1000 x 600	2	2	HC3206
160	1150 x 1000 x 600	2	4	HC3204
250	1550 x 1000 x 600	2	6	HC3203



ROWE CODE
HC3205

TOXIC SUBSTANCE

Indoor Storage Cabinets - **Metal**



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
30	770 x 515 x 465	1	1	HC3390
60	1065 x 515 x 465	1	2	HC3400
100	770 x 935 x 620	2	1	HC3009
160	1295 x 1115 x 500	2	2	HC0681
250	1825 x 1115 x 500	2	3	HC3305
250	1825 x 1440 x 585	2	3	HC0693



ROWE CODE
HC3009

OXIDIZING AGENT

Indoor Storage Cabinets - **Metal**



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
30	770 x 515 x 465	1	1	HC3403
60	1065 x 515 x 465	1	2	HC3405
100	770 x 935 x 620	2	1	HC0694
160	1295 x 1115 x 500	2	2	HC0695
250	1825 x 1115 x 500	2	3	HC0638



ROWE CODE
HC3405

PRATT
SAFETY SYSTEMS

SITE SAFETY SAFETY CABINETS

FLAMMABLE LIQUID

Indoor Storage Cabinets - Metal



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
30	770 x 515 x 465	1	1	HC3080
60	1065 x 515 x 465	1	2	HC3090
100	770 x 935 x 620	2	1	HC3100
160	1295 x 1115 x 500	2	2	HC3000
250	1825 x 1115 x 500	2	3	HC3200
350	1825 x 1440 x 585	2	3	HC3300
350	1825 x 865 x 865	1	3	HC0687
205	1825 x 865 x 865	1	Roller Set	HC0688
425	1825 x 1115 x 865	2	3	HC0689



ROWE CODE
HC3300

ORGANIC PEROXIDE

Indoor Storage Cabinets - Metal



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
30	770 x 515 x 465	1	1	HC0696
60	1065 x 515 x 465	1	2	HC0697
100	770 x 935 x 620	2	1	HC0698
100	1295 x 1115 x 500	2	2	HC0699
100	1825 x 1115 x 500	2	3	HC0700



ROWE CODE
HC0698

CLASS 4

Indoor Storage Cabinets - Metal



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
30	770 x 515 x 465	1	1	HC0701
60	1065 x 515 x 465	1	2	HC0702
100	770 x 935 x 620	2	1	HC0703
160	1295 x 1115 x 500	2	2	HC0671
250	1825 x 1115 x 500	2	3	HC0704



ROWE CODE
HC0702

DANGEROUS GOODS

Outdoor Storage Cabinets - Metal



VOL (L)	EXT DIMENSIONS (H x W x D mm)	DOORS	SHELVES	ROWE CODE
80	1375 x 785 x 635	1	2	HC0709
160	1375 x 785 x 900	1	2	HC0710
250	1375 x 1035 x 900	1	2	HC0711
250	1857 x 785 x 900	1	3	PC0530
350	1870 x 1185 x 900	1	3	HC0712
450	1870 x 1535 x 900	2	6	HC0669



ROWE CODE
HC0712

SITE SAFETY

GLOVE BOXES & DESICCATORS

Compact glove box, series 810-SG

PLAS ■ LABS, INC.™

With transfer chamber AND flat side access door.
Completely sealed glove box for basic isolation and/or containment.

The 810-SG Series by PLAS ■ LABS are constructed of an optically clear material, developed for isolation and / or containment. All units feature a high visibility of the interior, cleanliness, safety, complete containment, portability and a draft-free atmosphere. The installed purging valves are ideal for lowering oxygen and humidity levels.

OTHER FEATURES INCLUDE:

- Transfer chamber with purging valves and vacuum gauge
- Flat side access door for easily introducing larger equipment
- Four (4) key cock valves for purging:
 - Two on the main body, two on the transfer chamber
- Hospital grade multiple electrical outlet strip
- Pressure relief valve with small HEPA filter
- Standard 24 month warranty (not on gloves or consumables)
- White ambidextrous gloves.

Need to know more?
Need to know what other options
and models are available?
Need a quote? CALL US!



COMPACT GLOVE BOX

810-SG10

INSIDE DIMENSIONS 710W x 580D x 740Hmm

OUTSIDE DIMENSIONS 1100W x 610D x 790Hmm
(Inc transfer chamber 305mm long)

ACCESS DOOR OPENING 400W x 560Hmm

ROWE CODE

IG0182

***Other sizes available
upon request.***

Keep your samples SAFE with PLAS ■ LABS Desiccator Cabinets

DESICCATOR CABINET

ROWE
CODE

MODEL 860.5-CG - 300 x 300 x 150mm, clear acrylic with 2 shelves PD0086

MODEL 860-CG - 300 x 300 x 300mm, clear acrylic with 2 shelves PD0088

MODEL 860-CGA - 300 x 300 x 300mm, amber acrylic with 2 shelves PD0106

MODEL 861-CG - 300 x 300 x 600mm, clear acrylic with 4 shelves PD0084

MODEL 862-CG - 450 x 450 x 450mm, clear acrylic with 4 shelves PD0092

Plus many more to choose from



ROWE CODE
PD0106

SITE SAFETY

HAZARD WARNING LABELS & SIGNS

DANGEROUS GOODS CLASS LABEL FOR TRANSPORT	QTY	ROWE CODE
FLAMMABLE GAS CLASS 2.1 (flame symbol) - 20 x 20mm, self adhesive paper, red and white	Roll 250	HL0184
FLAMMABLE GAS CLASS 2.1 (flame symbol) - 50 x 50mm, red and white	250	HL1316
FLAMMABLE GAS CLASS 2.1 (flame symbol) - 100 x 100mm, self adhesive vinyl, black and red	25	HL0463
NON-FLAMMABLE NON-TOXIC GAS CLASS 2.2 (gas cylinder symbol) - 50 x 50mm, self adhesive paper, green and black	Roll 250	HL0167
NON-FLAMMABLE NON-TOXIC GAS CLASS 2.2 (gas cylinder symbol) - 100 x 100mm, self adhesive paper, green and black	250	HL0168
TOXIC GAS CLASS 2.3 (skull and crossbones symbol) - 100 x 100mm, self adhesive paper, black and white	Roll 1000	HL0474
FLAMMABLE LIQUID CLASS 3 (flame symbol) - 20 x 20mm, self adhesive paper, white and red	Roll 250	HL0418
FLAMMABLE LIQUID CLASS 3 (flame symbol) - 25 x 25mm, self adhesive paper, black and red	Roll 250	HL0187
FLAMMABLE LIQUID CLASS 3 (flame symbol) - 100 x 100mm, self adhesive vinyl, black and red	25	HL1203
FLAMMABLE SOLID CLASS 4.1 - 25 x 25mm, self adhesive paper, white, red and black	Roll 250	HL0162
OXIDISING AGENT CLASS 5.1 (flame over circle symbol) - 25 x 25mm, self adhesive paper, yellow and black	Roll 250	HL0185
TOXIC CLASS 6 (skull and crossbones symbol) - 25 x 25mm, self adhesive paper, black and white	Roll 250	HL0186
CORROSIVE CLASS 8 (corrosion symbol) - 20 x 20mm, self adhesive paper, black and white	Roll 250	HL1120
MISCELLANEOUS DANGEROUS GOODS CLASS 9 - 25 x 25mm, self adhesive paper, black and white	Roll 250	HL0453



ROWE CODE
HL0474



ROWE CODE
HL1203



ROWE CODE
HL1600

GHS HAZARD PICTOGRAM	QTY	ROWE CODE
EXPLOSIVE (exploding bomb symbol) - 15 x 15mm, GHS, black and red	Roll 100	HL0434
OXIDISING (flame over circle symbol) - 15 x 15mm, GHS, black and red	100	HL0435



ROWE CODE
HS0838

SAFETY / WARNING LABELS	QTY	ROWE CODE
CALIBRATION WRITE ON LABEL - 38 x 15mm, self adhesive vinyl cloth, green and white	350	HL1600
CALIBRATION BY DATE DUE WRITE ON LABEL - 38 x 15mm, self adhesive vinyl cloth, green and white	25	HL0195
CALIBRATED WRITE ON LABEL - 38 x 15mm, self adhesive vinyl cloth, black and white	350	HL0182
QUARANTINE WRITE ON LABEL - 70 x 100mm, removable self-adhesive, orange and black	Roll 500	HL0419
HAZCHEM - 100 x 100mm, self adhesive vinyl	250	HL0468
DANGER HAZARDOUS CHEMICALS - 100 x 100mm, self adhesive vinyl	250	HL0466
WARNING 1080 SODIUM FLUOROACETATE - 100 x 100mm, self adhesive vinyl	250	HL0467
SAFETY TAG - CAUTION, OUT OF SERVICE, DO NOT OPERATE - 125 x 75mm, non tear	100	HS0838



ROWE CODE
HS0819

SAFETY SIGNS	QTY	ROWE CODE
EMERGENCY SHOWER AND EYE WASH - 450 x 300mm, metal	Each	HS0819
HOT SURFACE - 250 x 180mm, self adhesive vinyl	Each	HS0613
DANGER HAZARDOUS CHEMICALS - 225 x 300mm, polypropylene	Each	HS0784
DANGER FLAMMABLE MATERIALS - 225 x 300mm, polypropylene	Each	HS0785
HAZCHEM - Three diamonds (Oxidizer 5.1 Corrosive 8 and Flammable liquid 3), 450W x 300mmH, metal, red on white with standard mounting holes	Each	HS0860



ROWE CODE
HS0613

SITE SAFETY

HAZARD CLEAN UP & SPILL STATIONS

HAZARD CLEAN UP - AUSTRALIAN MADE	SIZE	ROWE CODE
GLOBAL PEAT - O & F ABSORBENT Absorbs oil, fuel, coolants, degreaser, mild acids/bases, paint, ink and most chemicals Absorbent capacity: 30L per bag	28L Bag	NG0201
HAZCHEM ABSORBENT PADS - 30 X 30cm Absorbs acids, bases, coolants, paint, solvents, oil and hazardous chemicals Absorbent capacity: 16L per pack of 50	Pack 50	HH0055
HAZCHEM ECO TRUCK SPILL KIT Ideal for vehicles, forklifts, workshops, warehouses, factories and laboratories Absorbent capacity: Up to 35L		NS0027
HAZCHEM HIGH PERFORMANCE SPILL KIT Acids, bases, coolants, oils, fuels, solvents and hazardous chemicals Absorbent capacity: Up to 185L	120L	HS0796
HAZCHEM SPILL KIT Acids, bases, coolants, oils, fuels, solvents and hazardous chemicals Absorbent capacity: Up to 42L	20L Drum	HC0960
HYDROFLUORIC ACID SPILL KIT Absorbent capacity: Up to 650mL		HS0729
MERCURY DECONTAMINANT POWDER HgX	9kg	CM3772
MERCURY SPILL KIT- 3 x USE Absorbent capacity: components for 3 x spills	10L Pail	CM3921
OIL & FUEL ECONOMY PLUS SPILL KIT Suitable for unsealed areas in and around mine sites, drilling rigs and workshops Absorbent capacity: Up to 135L	120L	NS0026
OIL & FUEL ECONOMY PLUS SPILL KIT Oil, fuel, diesel, petrol, solvent and other petroleum products Absorbent capacity: Up to 235L	240L	HS0823
OIL & FUEL PADS - STANDARD DUTY - 45 X 45cm Use to absorb/clean up leaks, drips and spills or for wiping down tools, plant and machinery Absorbent capacity: 35L per pack	Pack 50	NA0051
SORBALITE ZEOLITE MINERAL SPONGE Absorbs oil, fuel, coolants, degreaser, mild acids/bases, paint, ink and most chemicals Absorbent capacity: 30L per bag	15kg	NS0028



ROWE CODE
NS0027



ROWE CODE
HS0796



ROWE CODE
CM3921



ROWE CODE
NS0028

HAZARD CLEAN UP	ROWE CODE
RELY ON VIRKON - HOSPITAL GRADE DISINFECTANT	
5g tablets (20 jars x 10 tablets, makes 500mL), TGA registered	CD1415
Shaker packs (6 x 500g, makes 50L), TGA registered	CD0591
5kg tub (makes 500L), TGA registered	CD1248
50g sachets (50 sachets, makes 5L), TGA registered	CD1249



ROWE CODE
CD0591

SITE SAFETY

HAZARD CLEAN UP & SPILL STATIONS

SPILL CONTAINMENT - AUSTRALIAN MADE

**ROWE
CODE**

FOUR DRUM SPILL BUND - POLYETHYLENE

Store up to 4 x 205L drums with removable PP grate to capture leaks, drips and spills. Load tested in Australia by NATA accredited engineers - dynamic compression test certificate available.

HP0251

IBC BUND - POLYETHYLENE

Store 1 x 1000L IBC to capture leaks, drips and spills. Forklift pockets provide easy relocation.

HI0027

TWIN IBC BUND - POLYETHYLENE

Store 2 x 1000L IBC's to capture leaks, drips and spills. Forklift pockets provide easy relocation. Store IBC's directly or palletised on the bund.

NT0054

BUNDED DECANT STAND - 40 LITRE

Use this bund to store and decant from 20L drums. Suits one or two 20L drums. 40L bund capacity.

HB0366



**ROWE CODE
HP0251**



**ROWE CODE
HI0027**



**ROWE CODE
NT0054**



**ROWE CODE
HB0366**



ROWESCIENTIFIC
PTY LTD www.rowe.com.au

If you cannot see what you are looking for, please call us

South Australia & NT
Ph: (08) 8186 0523

Queensland
Ph: (07) 3376 9411

Victoria & Tasmania
Ph: (03) 9701 7077

New South Wales
Ph: (02) 9603 1205

Western Australia
Ph: (08) 9302 1911

ORDER ONLINE AT YOUR CONVENIENCE

How to order online at your convenience using the Rowe Scientific Web Store and its easy to use features

As a guest, when you arrive at the Rowe Scientific web store, you can search for products and populate a trolley. To obtain pricing, please request a quote on your trolley of products. Clients with a login account will see pricing directly, and be able to place orders over the Rowe Scientific web store.

How to Log In

If you have a login account, enter your username and password in the login panel at the top right hand corner. This will log you into your personalised Web Store account and show your special pricing.

Your Welcome Panel

When you have successfully logged in, you will see your welcome panel to the right. Here you will find some useful features that are explained on this page.

My Backorders

My backorders is a searchable list of your current back orders with Rowe Scientific. This data is updated overnight to keep you fully informed.

Quick Order

YOUR TEN MOST PURCHASED ONLINE ITEMS			
Enter quantity on any of the items in the below lists then hit the 'Add all to Trolley' button at the bottom of the page.			
Items		Qty	UOM
CB3515	BUFFER pH 10.00 (BLUE) SOLUTION NIST TRACEABLE, IN CONTAMINATION PROOF BOTTLE (R608+)	\$19.71	500mL
FP1041	FILTER PAPER 25.0cm WET STRENGTHENED (MN722) IN BROWN PAPER PACKING	\$34.82	100

Quick Order allows you place orders from a list of routinely ordered products. There are two parts to your quick order page. Your top ten recent purchases or a manual list made up of items you add to your quick order list which are grouped by product type e.g. Glassware for glass beakers etc.

To add a product to your quick order list, click on the add to quicklist button on the product listing. If the button does not appear on an item, it is already in your quick order list.

ADD TO QUICKLIST

To order items directly by quick order list, simply enter the quantity required for each item being ordered then select the add to trolley button. You will see these items appear in your trolley panel and the total value of your order. You can continue searching and shopping for products or view your trolley and process your order.

Please note: only orders placed online appear in the list below	
view last 10 orders	
Order No. 211611 / Customer Ref. Caritest	
Product ID	Description
FF2227	FILTER PAPER 15.0cm SIMILAR TO Gr4 (MN617) (4340)
FF4648	FILTER PAPER GLASS FIBRE 11.0cm SIMILAR Gr GF/B

Order History

Here you can see a list of recent online orders placed via our web store.

The Search Bar

The search bar functionality on the right allows you to string key words together as you would a Google search. e.g. 'beaker 250mL' the results will list those items found which best match the exact criteria to the top of the results page. When you use the product filter buttons a second inner search bar appears in the search results page allowing you to drill down from the results using key words.

Search by YOUR part number

We can upload your part numbers (where supplied by you) for contract lines allowing you to also search by your part number for easy and accurate ordering of what you need. Your part number will appear in red in brackets directly after the Rowe code.

FF2900 (10617517) F
(MN1672) (472032)
Unit:100 Local Stock:64
Est. Lead days if no local or

See your agreed pricing

You will see your agreed contract and special pricing. This is indicated by the words **Your contract/special price** above the price.

Your Special price
\$10.00
VOLUME PRICING
QTY

Availability

You can check the availability and estimated lead times on products. Faster lead times may be available on request. The national stock is the total of our stock across the country.

TRACEABLE, IN CONTAMINATION PRO (R608+)		
Unit:500mL	Local Stock:463	National Stock:76
Est. Lead days if no local or national stock:5		
30S	C OF A	ADD TO QUICKLIST

www.rowe.com.au

South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au



www.rowe.com.au

The purpose of Rowe Scientific Pty. Ltd. is to keep our clients supplied with the scientific items they need to do their work.

HOW CAN WE HELP YOU TODAY?

NEW SOUTH WALES

Unit 1, 12 Essex Street
Minto, Sydney 2566
P (02) 9603 1205
E rowensw@rowe.com.au

QUEENSLAND

55 Westgate Street
Wacol, Brisbane 4076
P (07) 3376 9411
E roweqld@rowe.com.au

SOUTH AUSTRALIA

Lot 7, 41-47 O'Sullivan
Beach Rd, Lonsdale 5160
P (08) 8186 0523
E rowesa@rowe.com.au

VICTORIA

Unit 2, 42 Green Street
Doveton, Melbourne 3177
P (03) 9701 7077
E rowevic@rowe.com.au

WESTERN AUSTRALIA

11 Challenge Boulevard
Wangara, Perth 6065
P (08) 9302 1911
E rowewa@rowe.com.au