



ROWE SCIENTIFIC PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



END OF FINANCIAL YEAR SALE

Be Quick And Save On Marked Items - Offer Strictly Ends 30 June 2018 Or While Stocks Last

SYRINGE PUMP

TWO WAY PUSH PULL PROGRAMMABLE SYRINGE PUMP

- Supports two syringes simultaneously in a push pull situation
- Suitable for syringe size 10uL-60mL
- Flow rate 0.184uL/min – 83mL/min per syringe
- Linear speed 1um/min-150mm/min
- Linear force >20kgf
- High resolution 4.3" colour LCD, touch screen
- Five working modes
- Built-in RS485 communication & supports the MODBUS protocol
- Stainless steel case



ROWE CODE
IP0483

For further information and pricing contact your nearest Rowe Scientific representative.

MMM OVENS & INCUBATORS

VENTICELL® Oven Range

- The VENTICELL® oven range ensures a homogeneous temperature profile for all drying and heating processes due to its patented system of forced air circulation.
- 22, 55, 111, 222, 404 and 707 litre models (Internal volume of the chamber).
- Working temperature: 10 °C above ambient temperature up to 250 °C / 300 °C.



ROWE CODE
IO0128

5% OFF
ALL MMM OVENS

INCUCELL® Incubator Range

- The INCUCELL® line provides a soft air flow within the chamber. The fan forced variant features precise temperature distribution with minimal deviations.
- 22, 55, 111, 222, 404, and 707 litre models (Internal volume of the chamber).
- Working temperature: From 10 °C above ambient temperature up to 99.9 °C.



ROWE CODE
II0087

PARAFILM AND STOMACHER BAGS

Parafilm

In laboratories, Parafilm M® is a thermoplastic, self-sealing film that holds moisture loss to a minimum and offers excellent barrier protection to the contents of tubes, flasks, culture tubes and more.



ROWE CODE
PP1200
PP1250

GRADE Stomacher bags


These stomacher bags are primarily used in conjunction with blending machines to homogenise food samples.

The bags are consistently reliable, robust and have a negligibly low failure rate.

- GRADE Products are a world leading quality brand of laboratory stomacher bags.
- The bags are easy to open and use.
- Over 40 years of manufacturing know-how makes the GRADE range successful.



ROWE CODE
PB0119



ROWE CODE
PB0122



ROWE CODE
PB0121



ROWE CODE
PB0177

5% OFF
FROM THE ITEMS
LISTED BELOW

DRY ICE MAKER AND STORAGE

DILVAC dry ice maker

- Fast and efficient production of dry ice; 1kg blocks, on demand, in under 2 minutes in your own lab. Makes a 80mm x 110mm x 180mm block
- No electricity required
- Requires CO2 liquid gas cylinder with internal siphon
- Hose and gauge included as standard.

ROWE CODE
IDO103



DILVAC Dewar flasks

- Rust and corrosion resistant, enamelled outer container for durability and longevity
- Double walled and vacuum pumped borosilicate glass inner container
- for efficient storage, transport & use
- Carry handle and lid for ease of transportation
- Other models available

ROWE CODE
HD0052



WWW.ROWE.COM.AU

YOUR LAB + LIFE SCIENTIST PARTNER

YOUR 'ONE STOP SHOP' LABORATORY SUPPLY PARTNER

DLAB CENTRIFUGES

DLAB high speed laboratory centrifuges

DLAB D1008

- Small footprint
- For strips and 8x2ml tubes (2 rotors)
- 7000rpm motor

ONLY
\$299

ROWE CODE
IC0039



DLAB D3024

- Long lasting, high strength rotors - autoclavable
- High separation strength
- Max RCF up to 21380 xg
- 15000rpm motor

ONLY
\$1990

ROWE CODE
IC0036



BENCHTOP CENTRIFUGES

Benchtop refrigerated centrifuges

K243R refrigerated bench top centrifuge

- Large 3 litre capacity
- A range of rotors and adapters available including swing out and micro-plate rotors
- 500 - 15000rpm, capable of 22,000 G

ONLY
\$9900

ROWE CODE
IC4181



Consul 21 benchtop centrifuge

- 1.6 litre capacity
- A range of rotors and adapters for each size of tube
- 14300rpm
- 16 memory programs

ONLY
\$5200

ROWE CODE
IC0111



PIPETTE FILLERS

DLAB Levo manual pipette filler

- Designed for single handed operation, lightweight and easy to use.
- Precise pipetting control
- High applicability for most plastic and glass pipettes from 0.1mL to 100mL
- 3µm replaceable Hydrophobic filter

ROWE CODE
DP0099

ONLY
\$75



DLAB Levo Plus electronic pipette filler

- For most plastic and glass pipettes from 0.1mL to 100mL
- Designed for single handed operation
- Eight speeds of selection for aspiration and dispensing.
- Large LCD display showing low battery warning and speed setting
- Powerful rechargeable Li-ion battery enables 8 hours of runtime

ROWE CODE
DP0094

ONLY
\$235



BOECO & DLAB MICRO PIPETTERS

BOECO GP series premium adjustable volume pipette

- Completely autoclavable at 121°C
- Streamlined shaft for use in narrow reaction vessels
- Colour-coded display window for easy tip selection



| ROWE CODE | VOLUME |
|-----------|----------------|
| DP4924 | 0,5 - 10 µl |
| DP2117 | 2-20 µL |
| DP2121 | 10-100 µL |
| DP2122 | 20-200 µL |
| DP2131 | 100-1000 µL |
| DP2132 | 500-5000 µL |
| DP2166 | 1000-10,000 µL |

DLAB mechanical adjustable volume pipette

- Fully autoclavable
- Clear digital display
- Supplied with individual calibration certificate



 DLAB MECHANICAL PIPETTERS
\$155

| ROWE CODE | VOLUME |
|-----------|-----------|
| DP0114 | 2 - 20 µl |
| DP0116 | 10-100 µL |
| DP0104 | 20-200 µL |
| DP0119 | 1-5 mL |

Full range available.

POLYSCIENCE BENCHTOP CHILLERS

LS, LM and MM series of Polyscience chillers

These powerful, low-temperature chillers are well matched for use with rotary evaporators, vacuum systems, spectrometers, and other analytical equipment. They also are available with a mobile cart accessory for convenient placement under a bench.

Features

- Large, easy to read LED display
- Space-saving design
- Cooling at ambient temperatures as high as 35°C
- Low flow shut-off and alarm, high and low temperature alarms
- Front mounted fluid level gauge
- Simple setup, operation, and maintenance
- Choice of pumps



ROWE CODE
IC0273

ROWE CODE
IC0275

ROWE CODE
IC0270

Other models available on request, get in touch with us now.

VELP BOD SYSTEMS

VELP BOD Sensors

- System 6 and system 10
- Available with a wireless databox system
- Option to purchase individual BOD sensor with amber bottle
- Many options available to suit your requirements



ROWE CODE
IB6050



ROWE CODE
NV4100

VELP Cooled incubator 120 & 215 litres

- Accepts BOD sensor system 6, or System 10
- Internal sockets to power BOD measuring devices
- **II1643** comes with selectable temperature ranges from 3 to 50°C

ROWE CODE
IV4004

ROWE CODE
II1643



Ask us to recommend a system for your laboratory.

TARSONS WIDE MOUTH REAGENT BOTTLES

- Tarsons wide mouth reagent bottles available in polypropylene and HDPE
- Complete with screw cap which are leak-proof even without the aid of a cap or gasket.
- Ideal for liquid storage and shipping liquid samples
- 6 sizes available

ORDER OVER \$200 FROM THE LIST BELOW AND TAKE **10% OFF**



| ROWE CODE | VOLUME | MATERIAL TYPE |
|-----------|--------|---------------|
| PB0204 | 30 mL | PP |
| PB0207 | 60 mL | PP |
| PB3335 | 125 mL | PP |
| PB3337 | 250 mL | PP |
| PB3319 | 500 mL | PP |
| PB3390 | 1 L | PP |

PP shown, HDPE available

Also available, Tarsons beakers, filter holders, vials and racks, funnels, jugs and carboys with stopcocks.

SIMAX LAB BOTTLES

SIMAX laboratory bottles

Reagent bottles with screw top cap,

Features include

- Borosilicate glass
- Red and blue caps available
- pouring rings to suit
- marked graduations

SIMAX



10% OFF
ALL SIMAX GLASSWARE

| ROWE CODE | VOLUME |
|-----------|---------|
| GB7910 | 100 mL |
| GB7955 | 250 mL |
| GB8005 | 500 mL |
| GB8055 | 1000 mL |
| GB8105 | 2000 mL |

FILTER-BIO SYRINGE FILTERS

Filter-bio syringe filters

Glass fibre pre filters available

FILTER-BIO



Membrane types including

- PES (polyethersulfone)
- PVDF (polyvinylidene difluoride)
- MCE (mixed cellulose esters)
- CA (cellulose acetate)
- CN (cellulose nitrate)
- Nylon, polypropylene
- PTFE (polytetrafluoroethylene)
- polycarbonate and glass-fibre prefilter.
- Membranes are available in a huge range of diameters and pore sizes.
- A large variety of samples are available for trial.



ROWE CODE **FF0234**

- 25 mm
- 0.45 micron
- Pack of 100
- Polyethersulfone



ROWE CODE **FF0094**

- 30 mm
- 0.45 micron
- Pack of 1000
- Polyethersulfone

SIMAX GLASSWARE

Rowe Scientific stocks a range of SIMAX glassware in beakers, centrifuge tubes, desicators and flasks

- made from borosilicate glass
- smooth and non porous, perfectly transparent, indifferent to catalysts, corrosion resistant even in a heavy duty operation of up to 300 degrees Centigrade

10% OFF
ALL SIMAX GLASSWARE



ROWE CODE **GD1219**



ROWE CODE **GF1855**



ROWE CODE **GB4085**

TARSONS CRYOGENIC VIALS

Polypropylene cryogenic vials with screw caps

Cryovial



ROWE CODE **PV1030**

ROWE CODE **PV1025**

| ROWE CODE | VOLUME |
|-----------|--------|
| PV1025 | 1.8 mL |
| PV1030 | 4.5 mL |

Cryochil freezer boxes

ROWE CODE **PT2005**

- Holds 0.5, 1.5, 2.0 mL tubes, 1.0 and 1.8 mL cryo vials
- autoclavable
- hinged lid



BONUS OFFER

Buy 10 packs of the cryogenic vials, and receive a free PT2005 Cryochill freezer box rack valued at \$16

TARSONS FILTER TIPS

MAXIPENSE™ Ultra Low Retention filter tips

Why Tarsons tips?

- Medical grade resin conforming to USP Class VI. Sterile and ready to use
- Made using 100% premium virgin polypropylene, free from heavy metals, lubricants and dyes - Autoclavable
- Manufactured in a fully automated class 100K clean room making Tarsons tips certified DNase, RNase and pyrogen free
- Tips contain a hydrophobic polyethylene filter which protects samples and pipettes from contamination
- Especially suited to sensitive PCR and real time PCR applications.



5% OFF
FROM THE ITEMS PURCHASED BELOW

| ROWE CODE | TIP CAPACITY | LENGTH |
|-----------|--------------|--------|
| DT0025 | 10µL | 30mm |
| DT0026 | 20µL | 53mm |
| DT0027 | 50µL | 55mm |
| DT0028 | 100µL | 55mm |
| DT0029 | 200µL | 55mm |
| DT0030 | 300µL | 55mm |
| DT0031 | 1000µL | 76mm |

Pack size:

10 racks x 96 tips (all graduated)

Brand Compatibility:

Includes Eppendorf®, Biohit® (Boeco), Finnpiptette®, Gilson®, Rainin®, Nishiryō®

HPS AND MERCK STANDARDS

Rowe Scientific supplies a comprehensive range of element standards from HPS and Merck. These lines offer the choice of mineral acid matrix, providing options to support a wide range of methods and applications, including adherence to the EU Pharmacopoeia monographs.

High Purity Standards
Accredited to ISO17025 and ISO Guide 34

- IC Standards
- ICP and AAS
- Metallo Organic
- Pharmaceutical standards
- Single element, Multi's and Custom
- Organic reference standards



Contact your nearest Rowe Scientific representative for more information.

CHEMICALS AND STANDARDS

Rowe Scientific supply a range of reagents, standards and fine chemicals to accommodate your laboratory requirements.



ULTRA Scientific Chemicals and standards

- 1mL ampoules to 20L drums
- Pesticides
- PCB's plus related compounds
- Lab and method check standards



Contact your nearest Rowe Scientific representative or visit www.rowe.com.au

CITOTEST MICROSCOPE SLIDES

CITOTEST® have a range of slides for your specific laboratory microscope application.

Citotest slides come in

- Frosted
- ColorCoat™
- ColorInx™
- PrintAID™
- ColorPlus™
- AutoSmear™
- PapSmear™
- Latex glass
- Slides with cavity
- Cover glass
- BoroGlas™
- EstaGlas™
- Circular Cover Glass
- Hemacytometer Cover Glass



ROWE CODE
GM0006



5% OFF
ALL CITOTEST
MICROSCOPE
SLIDES

Contact your nearest Rowe Scientific representative today to take advantage of this great offer.

SEROLOGICAL PIPETTES - STERILE



Buy 1000 serological pipettes and receive a bonus HR0021 storage rack for free.



1 per Customer.

BONUS OFFER

Bonus HR0021 storage rack with orders over 1000



ROWE CODE
HR0021

| ROWE CODE | CAPACITY | UOM | PRICE |
|-----------|----------|-----|-------|
|-----------|----------|-----|-------|

Individually Wrapped

| | | | |
|--------|-------|------|-------|
| PP0001 | 1 mL | 1000 | \$185 |
| PP0002 | 2 mL | 800 | \$150 |
| PP0003 | 5 mL | 300 | \$99 |
| PP0004 | 10 mL | 200 | \$84 |
| PP0005 | 25 mL | 200 | \$130 |
| PP0006 | 50 mL | 100 | \$109 |

Bulk Packed

| | | | |
|--------|-------|------|-------|
| PP0019 | 1 mL | 1000 | \$153 |
| PP0020 | 2 mL | 800 | \$143 |
| PP0021 | 5 mL | 500 | \$125 |
| PP0022 | 10 mL | 350 | \$100 |

GLOVES

Rowe Scientific has a range of nitrile, neoprene and powder free latex gloves to suit many different applications within the laboratory.

Choose from Ansell, Supermax or RCR for your next bulk glove purchase.

Ansell Touch N Tuff



ROWE CODE
HG1893

Supermax powder free nitrile



ROWE CODE
HG1908

RCR powder free



ROWE CODE
HG1838

For information on any of these products and many more, please contact your nearest Rowe Scientific consultant, or visit our online webstore now.



ROWE SCIENTIFIC
PTY LTD www.rowe.com.au
For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987. We are proudly a 100% Australian owned company.



South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

New South Wales & ACT
Ph: (02) 9603 1205
rowesw@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

www.rowe.com.au