



MAXIPENSE™ Ultra Low Retention

Filter Tips



EXAMPLES **INSIDE**



ROWE SCIENTIFIC

PTY LTD www.rowe.com.au

Please find a selection of Filter tips provided to identify the best fit for your pipette. Take note of the **Rowe Code** to reorder from your Account Manager or order online at **www.rowe.com.au**

Sterile racked samples of MAXIPENSE™ tips are available to evaluate as required.

Why Tarsons Tips?

1. Medical grade resin conforming to USP Class VI. Sterile and ready to use
2. Made using 100% premium virgin polypropylene, free from heavy metals, lubricants and dyes - Autoclavable
3. Manufactured in a fully automated class 100K clean room making Tarsons tips certified DNase, RNase and pyrogen free
4. Tips contain a hydrophobic polyethylene filter which protects samples and pipettes from contamination
5. Especially suited to sensitive PCR and real time PCR applications.

MAXIPENSE™ filter tips are supplied in racks-sterile.

Standard filter tips are available either in bulk bags or racked - sterile.

Please request samples from your local Rowe Scientific representative or office in your state.



Product Specifications:



MAXIPENSE™

Low retention surface is achieved by treating the polypropylene at the molecular level, modifying the surface of the polymer to have a very low co-efficient of friction.

Pack size: 10 racks x 96 tips (all graduated)
Brand Compatibility: Includes Eppendorf®, Biohit® (Boeco), Finnpipette®, Gilson®, Rainin®, Nishiryo®

Rowe Code: DT0025
Tarsons Code: T529100
Tip Capacity: 10µL
Length: 30mm



Rowe Code: DT0026
Tarsons Code: T529101
Tip Capacity: 20µL
Length: 53mm



Rowe Code: DT0027
Tarsons Code: T529102
Tip Capacity: 50µL
Length: 55mm



Rowe Code: DT0028
Tarsons Code: T529103
Tip Capacity: 100µL
Length: 55mm



Rowe Code: DT0029
Tarsons Code: T529104
Tip Capacity: 200µL
Length: 55mm



Rowe Code: DT0030
Tarsons Code: T529105
Tip Capacity: 300µL
Length: 55mm



Rowe Code: DT0031
Tarsons Code: T529106
Tip Capacity: 1000µL
Length: 76mm





Offering you competitively priced precision pipettes from the World's Leading Scientific Brands.



Contact us today about solutions to your pipetting needs.



www.rowe.com.au



Quality Endorsed Company
 ISO 9001:2008
 UK 10372
 SA 06304

South Australia & NT
 Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
 Ph (07) 3376 9411
rov/eqld@rowe.com.au

Western Australia
 Ph: (08) 9302 1911
rowewa@rowe.com.au

New South Wales
 Ph: (02) 9603 1205
rowensw@rowe.com.au

Victoria & Tasmania
 Ph: (03) 9701 7077
rowevic@rowe.com.au