



ROWE SCIENTIFIC PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



MAKE THE DISCOVERY.

EXPLORE BEYOND THE LABORATORY OVEN



STERICELL®

HOT AIR STERILISER

SPECIAL OFFER

ONLY \$3499

While Current Stocks Last.

IS0174

Volume: **111L**

Temp: **Ambient +10°C - 250°C**

Interior: **Stainless Steel**



STERICELL® is ideal for hot air sterilisation of materials and instruments whilst maintaining a stable temperature for a pre-determined duration.

The built in fan quietly circulates air within the chamber to ensure an even distribution of heat and eliminates the formation of "cold air" areas within the chamber.



We have product specialists ready to help you choose your next oven

www.rowe.com.au

ECOCELL®

NATURAL AIR CIRCULATION OVEN



A line of economical ovens with a wide temperature range, which reliably heats and dries materials to exact temperatures.

The ECOCELL® line produces no noise and provides a very soft air convection within the chamber to uniformly heat materials within the oven.

I00134 55L **I00133** 111L **I00132** 222L

FRIOCELL®

FORCED AIR INCUBATOR WITH ACTIVE COOLING



The FRIOCELL® incubator allows precise tempering of samples.

The unique cooling system offers an exact simulation of a natural process while reducing the samples evaporation rate.

I10059 55L **I10060** 111L

VACUCELL®

VACUUM CHAMBERED OVEN



Temperature sensitive and oxidative materials such as powders and granulated samples can be safely dried in the VACUCELL® vacuum drying oven. This is achieved by either evacuating air with an external vacuum connection or by injecting inert gas to displace the air within the chamber.

This oven is also ideal for drying intricate components and hardware with inaccessible holes and deep crevices.

I00111 55L

South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au