



pH CUBE



The pH Cube meter is a compact benchtop unit which is perfect for education and QA labs.

The pH Cube provides accurate and reliable measurements of pH, ORP and temperature, and are simple and very easy to use.



Tough, all steel construction
Bright LED display
Manual controls
Easy to use
Choice of 3 colours
(Red, Blue and Silver)



MODE	RANGES	RESOLUTION	ACCURACY
Temp	-10.0 to 120.0°C	0.1°C	+/-0.2°C
pH	0.00 to 14.00	0.01 pH	+/- 0.01 pH
ORP (mV)	0 to +/-1999mV	1 mV	+/-1 mV
Input Impedance	> 3 x 10 ¹² 10 Ω		
Asymmetry Range	-1.00 to +1.00 pH		
Slope Range	85 to 105%		
Temperature Compensation	0 to 100°C Automatic or Manual		
Environment	0 to 45°C temperature 0-90% R.H.		
Dimensions	123 x 123 x 123mm		
Weight	850g		
Power	12V DC 200mA via AC/DC adaptor.		



South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

www.rowe.com.au