

# get the **burkle** SCOOP!



www.rowe.com.au



South Australia &  
Northern Territory  
Ph: (08) 8186 0523  
rowesa@rowe.com.au

Queensland  
Ph: (07) 3376 9411  
roweqld@rowe.com.au

Victoria & Tasmania  
Ph: (03) 9701 7077  
rowevic@rowe.com.au

Western Australia  
Ph: (08) 9302 1911  
rowewa@rowe.com.au

New South Wales & ACT  
Ph: (02) 9603 1205  
rowensw@rowe.com.au

SCOOPS  
with or without lids  
Individually packed  
25, 50, 100, or 150mL




SPOONS  
Individually packed  
2.5 or 10mL



**NEW**  
FROM ROWE  
SCIENTIFIC

SPATULA  
with or without cover  
Individually packed



Quality  
Endorsed  
Company  
ISO 9001:2008  
LIC 10372  
SAI Global

- Individually packed
  - Gamma sterilized
  - Range of sizes
- Ideal for the pharmaceutical, food & feed laboratories and more.

**SPECIAL BULK BUY PRICING AVAILABLE**



The **LaboPlast**<sup>®</sup> manufactured under clean room standards and **SteriPlast**<sup>®</sup> Gamma sterile range of single-use scoops, spoons and spatulas are made from polystyrene and produced, assembled and individually packed in accordance to clean-room class 7 standards, and in within EU & FDA guidelines.

## LaboPlast<sup>®</sup>

## SteriPlast<sup>®</sup> Gamma sterile

ROWE CODE	SCOOPS	PCS PER PACK	LIST PRICE PER PACK	BUY 10 PACKS AND OVER PER PACK	ROWE CODE	PCS PER PACK	LIST PRICE PER PACK	BUY 10 PACKS AND OVER PER PACK
PS0063	SCOOP SINGLE USE WITHOUT LID 141mm 25mL PLASTIC	10	30.67	27.56	PS0048	10	33.34	29.96
PS0065	SCOOP SINGLE USE WITH LID 142mm 25mL PLASTIC	10	47.94	43.08	PS0049	10	50.73	45.65
PS0066	SCOOP SINGLE USE WITHOUT LID 160mm 50mL PLASTIC	10	31.63	28.43	PS0051	10	34.43	32.36
PS0067	SCOOP SINGLE USE WITH LID 161mm 50mL PLASTIC	10	49.39	44.39	PS0052	10	52.19	46.90
PS0068	SCOOP SINGLE USE WITHOUT LID 205mm 100mL PLASTIC	10	33.58	30.18	PS0053	10	35.64	32.04
PS0069	SCOOP SINGLE USE WITH LID 206mm 100mL PLASTIC	10	51.95	46.69	PS0054	10	53.66	48.22
PS0070	SCOOP SINGLE USE WITHOUT LID 216mm 150mL PLASTIC	10	33.58	31.55	PS0055	10	36.48	32.79
PS0071	SCOOP SINGLE USE WITH LID 218mm 150mL PLASTIC	10	54.30	48.81	PS0056	10	54.87	49.32

The **LaboPlast**<sup>®</sup> manufactured under clean room standards and **SteriPlast**<sup>®</sup> Gamma sterile range of scoops are available in four sizes. They're also available with transparent cover to secure a sterile and contamination-free sample. The handle construction and slightly angled edge allows for the scoop to be laid down without the material trickling out. The scoops are ideal for sampling powder, granules and pastes.




ROWE CODE	SPOONS	PCS PER PACK	LIST PRICE PER PACK	BUY 10 PACKS AND OVER PER PACK	ROWE CODE	PCS PER PACK	LIST PRICE PER PACK	BUY 10 PACKS AND OVER PER PACK
PS0072	SPOON SINGLE USE PACKED 127mm 2.5mL PLASTIC	100	73.79	66.37	PS0057	100	77.06	69.12
PS0073	SPOON SINGLE USE PACKED 170mm 10mL PLASTIC	100	78.72	70.74	PS0058	100	82.20	73.73

The **LaboPlast**<sup>®</sup> manufactured under clean room standards and **SteriPlast**<sup>®</sup> Gamma sterile range of spoons are made of polystyrene are produced in a clean-room and are individually packed. The sampling spoon is available in 2.5 ml (similar to teaspoon) and 10 ml (similar to tablespoon). Ideal for sampling of powder, granules and pastes.



ROWE CODE	SPATULAS	PCS PER PACK	LIST PRICE PER PACK	BUY 10 PACKS AND OVER PER PACK	ROWE CODE	PCS PER PACK	LIST PRICE PER PACK	BUY 10 PACKS AND OVER PER PACK
PS0074	SPATULA SINGLE USE W/O COVER 263mm INSERT TO 150mm PLASTIC	10	29.44	26.45	PS0059	10	36.33	27.57
PS0075	SPATULA SINGLE USE WITH COVER 263mm INSERT TO 150mm PLASTIC	10	44.29	39.80	PS0061	10	46.25	41.48



The **LaboPlast**<sup>®</sup> manufactured under clean room standards and **SteriPlast**<sup>®</sup> Gamma sterile range of Spatulas feature a long ergonomic solid handle & pointed blade means it can be easily pierced into paper and plastic bags. Ideal for sampling of powder, granules & pastes. As an option the spatula is available with cover in order to secure a sterile and contamination free sampling process.



Also available from the **burkle** sampling range ...



### Measuring scoops. Ref:96140002

Universal graduated scoops ideal for measuring powders and granules. For exact metering, overfill then wipe the excess at the upper edge.


- Available from 2mL to 1000mL capacity
- Polypropylene, white
- High chemical-resistance
- Calibrated volume when full to the brim



### Pharma/Food scoops. Ref:53240010

Like all samplers for the pharmaceutical industry, they are made of Stainless steel V4A without edges, crevices or undercuts in which sample residues could collect and contaminate subsequent samples.

The highly polished inner and outer surfaces of the scoop meet all the requirements for GMP compliance in production and sampling.



### Zone samplers. Ref:53172055

The big advantage of zone samplers is that all-layer samples can be collected from all layers of the container. Suitable for various bulk goods from the finest powder to large grains such as corn or nuts.



Samplers are available in different lengths.

55 cm suitable for sacks & bulk bags. 85 cm (for barrels, tubs, mixed containers)

150 cm and 250 cm (Jumbo sampler) for big bags, silos, and tanks.

Chamber combinations illustration above: **1.** Multi (closed chambers), **2.** Uno (one closed chamber) and **3.** All-layer sampler (open inner tube)



### Pharma/Food Close it labels. Ref:53033017

Burkle have developed close-it food, a control seal with an adhesive which complies with the European guidelines for foodstuffs. \*Close-it food reliably prevents products from leaking out after sampling.

It has a barrier layer of aluminium foil which makes it suitable for closure after sampling of foodstuffs, pharmaceuticals, cosmetics, food additives, even frozen goods. In short, for all applications in which cleanliness is particularly important.

Available 95mm square or 150mm square in red, yellow, blue, green black and white. Particularly suited for the food and pharma industry.

\*FDA approval certificate available for all products