



A 100% Australian owned company supplying scientific laboratories throughout Australia and overseas since 1987.



SEPTEMBER 2021

Discover some of the latest and most popular laboratory products available from Rowe Scientific Pty Ltd.

July

July 1 - 1934, the first X-ray photograph of the whole body taken by Arthur W. Fuchs of the Eastman Kodak Company.



August

National Science Week
August 14 - 22

September

Sept 25 - 1974, scientists first reported that freon gases released from aerosol spray cans were destroying the ozone layer.



Source: Today In Science (www.todayinsci.com). Images Shutterstock.

KN95 Face masks

With ear loops.
Splash resistant Type IIR

This disposable white KN95 face mask is soft and offers a comfortable and easy to breathe through user experience.

- Made from non-woven, static cotton, melt-blown fabric
- Individually packed
- Professionally tested in Australia
- 3 layers of protection
- 99.9% bacterial filtration efficiency
- TGA listed (ARTG #336093)



HM0212
25 per pack

\$25
Per pack of 25



- ✓ One stop shop
- ✓ 100% Australian
- ✓ Online ordering
- ✓ Phone ordering
- ✓ Account Managers
- ✓ Leading brands
- ✓ Chemicals
- ✓ Consumables
- ✓ Equipment



Quality Endorsed Company
ISO 9001:2015
SAI Global



Bel-Art Products



Gilian

Haier

Kartell LABWARE



www.rowe.com.au

Quantumclean Concentrate

For total manual cleaning (including radioactive decontamination) of all laboratory and industrial glassware, plastics, ceramics and stainless steel.
Australian made by Rowe Scientific Pty Ltd.



CQ1000

CQ1050

CQ1100

PART PROCEEDS DONATED TO THE ROYAL FLYING DOCTOR SERVICE OF AUSTRALIA. FIND OUT MORE ON PAGE 8

VELP FLOCCULATORS

Do you test wastewater for coagulation and turbidity?
Designed for use with the Jar test, a preliminary test used in
water treatment and similar applications.



IF1160 FP4 Portable Flocculator

4 Place Portable Jar Test
Flocculator, 240v AC 12v DC

- Choice of 20, 40, 50, 100 and 200 rpm stirring speeds
- Precise speed from an electronically controlled DC gear motor
- 4 samples can be processed in parallel at one time
- Designed to be used on-site for immediate results
- Strong resistance to chemical and mechanical corrosion
- Central light to facilitate reading.



Also available,
sold separately
IF1165
Flocculator
Carrying Case
for FP4



IF1171 FC S Series Flocculator with 4 or 6 independent positions

Flocculators with 4 or 6
independent positions that
allow standard conditions to be
adopted for the Jar Test.

- Speed can be changed whilst in operation
- Illuminated back panel to aid in sample observation
- Stirring rods adjustable in height without tools
- Precise speed from an electronically controlled DC gear motor
- Strong resistance to chemical and mechanical corrosion
- Independently selectable speeds 10, 15, 30, 45, 60, 90, 120, 150, 200, 300 rpm.



IF1155 JLT 4 Series Flocculator with 4 or 6 positions

Digital Flocculators with 4
and 6 positions with timer
and illuminated back panel
to allow standard conditions
to be met for the Jar Test.

- Microprocessor controlled operating time for unattended operation
- Illuminated back panel to aid in sample observation
- Stirring rods adjustable in height without tools
- Digitally selectable speed from 10 to 300 rpm
- Identical speed and time for all positions for comparable results
- Strong resistance to chemical and mechanical corrosion
- High visibility LED displays for speed and time.



A 100% Australian owned
company since 1987.

ROWE SCIENTIFIC PTY LTD

LABLIKES

ISSUE 11

VELP OVERHEAD STIRRERS



Ermes enabled

The ADVANCE models come with a 3.5" colour display which enables quick method setting and displays precise torque and temperature values. Equipped with a vibration sensor, timer function and Wi-Fi to enable connection to the VELP Ermes cloud platform for remote operation and monitoring.

OHS ADVANCE

SAFETY

IP54 rating for protection against corrosion and harsh environments

The "Lock key" prevents unintentional setting changes

PERFORMANCE

SpeedServo™ technology with a high powered brushless motor ensures constant torque at any speed

Your choice of single or double speed models

Up to 200 Ncm torque

HIGH TECH

The OHS SmartChuck™ technology enables a one handed and tool free changing of shafts

Included is a PT100 temperature probe

OHS DIGITAL

SAFETY

IP54 rating for protection against corrosion and harsh environments

The "Lock key" prevents unintentional setting changes

Over current and over temperature protection

PERFORMANCE

SpeedServo™ technology with a high powered brushless motor ensures constant torque at any speed

Your choice of single or double speed models

Up to 200 Ncm torque and 100 L volume depending on the model

'Silent' operation even at maximum speed (up to 2000 rpm)

HIGH TECH

The OHS SmartChuck™ technology enables a one handed and tool free changing of shafts

Digital display for precise setting and monitoring

Integrated timer for unattended operation

Real time torque LED indicator

USB interface
Included is a PT100 temperature probe

Beaker Borosilicate Glass

BRANDS WE CARRY

Simex, Kimble heavy duty, BRAND, Witeg volumetric. Please enquire for more.

SIMAX



Kimble

witeg

	mL	PACK SIZE	ROWE CODE
Low Form Beakers	10	ea 10	GB1509
	25	ea 10	GB1560
	50	ea 10	GB1660
	100	ea 10	GB1810
	150	ea 10	GB1870
	250	ea 10	GB2060
	400	ea 10	GB2210
	600	ea 10	GB2310
	800	ea 10	GB2510
	1000	ea 10	GB2610
	2000	ea 4	GB2810
	3000	ea 4	GB2960
	5000	ea 2	GB3110
	10000	ea	GB0356
	20000	ea	GB0357
Heavy Duty	250	ea 48	GB2100 GB2102
	400	ea 48	GB2340 GB2342
	600	ea 36	GB2350 GB2351

	mL	AVL or PACK SIZE	ROWE CODE
	25	ea 10	GB0358
	50	ea 10	GB1366
Tall Form Beakers	100	ea 10	GB3210
	150	ea 10	GB3310
	250	ea 10	GB3460
	400	ea 10	GB3560
	600	ea 10	GB3850
	1000	ea 10	GB4085
	2000	ea 6	GB4071
	3000	ea 6	GB4065 GB0359



STIRRING PADDLES AVAILABLE AT ADDITIONAL COST



6-hole paddle



Anchor



Fixed blade



Folding blade



Floating blades



Propeller



Turbo propeller



Turbine

All special pricing expires 31 October 2021 or until stocks last and is subject to GST.

Purchase the New Brand® HandyStep® touch S

at \$2,190.00

from ROWE SCIENTIFIC (Rowe Code: DH0031)

AND

receive the
charging stand valued at
\$654.36, **Free** of charge!

(Rowe Code: DC0046)

**FREE
CHARGING
STAND**

OFFER EXPIRES
30 SEPTEMBER
2021



Versatile

HandyStep® touch with the standard modes Multi-Dispensing, Auto-Dispensing, Pipetting and Favorites;

HandyStep® touch S with the additional functions Sequential Dispensing, Multi-Aspiration and Titration

Long-lasting

Many hours of continuous operation with a single battery charge

Robust

High-quality materials provide long service life, and smooth surfaces make cleaning / disinfection easy

Efficient

In addition to BRAND PD-Tips, compatible third party tips can also be used

Easy

Automatic size detection of BRAND PD-Tips II

Simple

Automatic tip ejection by the press of a button prevents contamination

USB-C connection

For quick charging by cable

STEP button

Ergonomically arranged STEP button for effortless volume filling and delivery





A 100% Australian owned company since 1987.

ROWE SCIENTIFIC PTY LTD

LABLIKES

ISSUE 11



Easy to understand

Main menu with all modes for quick selection



Clearly structured

Working mode with all important information



Organized

All options visible at a glance



Helpful

Internal help right where it is needed

The perfect team

Automatic size detection with PD-Tips II from BRAND

The HandyStep® touch saves time and prevents errors through automatic tip size recognition of the PD-Tips II from BRAND, with patented coding on the pistons. After inserting the tip, the size is automatically recognized and displayed, making it easy to select the volume to be dispensed. When a new PD-Tip II of the same size is inserted, all instrument settings are maintained.



Multi dispensing with the HandyStep® touch and PD-Tip II



Sequential dispensing with the HandyStep® touch S and PD-Tip II

VOLUME mL	Pack size	NON STERILE ROWE CODE	BIO-CERT ROWE CODE
0.1	100	CP3095	DH0012
0.5	100	DH0020	DH0013
1	100	DH0021	DH0014
1.25	100	CO3096	DH0015
2.5	100	DL3430	DH0016
5	100	DL3445	CP3104
10	100	DH0025	DH0017
12.5	100	PT0344	CP3097
25	50+	DH0024	DH0018
50	25+	CP3102	DH0019

PD TIP Selection guide



PROTECTIVE EYEWEAR

bolle Quality safety glasses for less than \$5 a pair (20 pair min)

BL30 B LINE

B-LINE is a range of safety eyewear adapted to common industrial risks within the laboratory and beyond.

With its ultra wrap-around fit, the BL30 range meets the requirements for many work environments, including laboratory chemists, first responders and frontline health workers. Equipped with an anti-scratch and anti-fog coating on both sides of the lens, BL30 offers long-lasting comfort.

- ➕ ANTI-FOG & ANTI-SCRATCH COATING ON BOTH SIDES OF THE LENS
- ➕ LIGHT WEIGHT FRAME FOR ALL DAY COMFORT
- ➕ BUILT IN STOP-HINGE SYSTEM (ON THE HINGE OF THE TEMPLES)
- ➕ UNIVERSAL FIT
- ➕ POLYCARBONATE HINGE WITH NO SCREWS FOR ADDED RELIABILITY

bolle
SAFETY



All special pricing expires 31 October 2021 or until stocks last and is subject to GST.

Miele

Greater capacity. Greater purity. Greater flexibility.

Laboratory glassware reprocessing for analytical purposes with the PG 85 series

Miele Professional. Immer Besser.

Purchase any of the below Miele professional glassware washers and receive **FREE*** glassware

• Order a PG8504 and receive \$500.00 worth of glassware **FREE***

• Order a PG8583AD or PG8583ADP and receive \$1000.00 worth of glassware **FREE***

• Order a PG8583 CD ADP or PG8583ADP and receive \$1500.00 worth of glassware **FREE***

*FREE glassware must be selected from a supplied list and pack size, only while stocks last. Only available on glassware washer models listed. Contact your local Rowe Scientific consultant to find out more.



Model number PG8504

Rowe Code **CD1543**

PG8583 AD

ID0236

PG8583 ADP

ID0104

PG8583 CD ADP

IG0046

PG8593 AD

IG0047



A 100% Australian-owned
company since 1987.

ROWE SCIENTIFIC PTY LTD

LABLIKES

ISSUE 11

TUBE RACKS

AUTOCLAVABLE POLYPROPYLENE

PR1117

18mm diameter 40 place (4x10) 36mm hole depth **low-profile** solid red glass filled polypropylene temp range (freezer -20 autoclave +150c)




PR0171

12mm hole x 60 place polypropylene autoclavable red




PR0172

12mm hole x 60 place polypropylene autoclavable green



PT0448

18mm diameter 40 place (4x10) 36mm hole depth low-profile solid green glass filled polypropylene temp range (freezer -20 autoclave +150c)




PR1115

18mm 40 place (4x10) solid red glass filled polypropylene temp range (freezer -20c autoclave +150c)




PR1114

13 and 14mm tubes 60 place (5x12) solid red polypropylene



PR0025

16mm hole 44 place (4x11) numbered temp range (freezer -20c autoclave +150c) solid red polypropylene



Label & tape solutions for today's laboratories



HL0278 \$741.00
Printer labels for Liquid Nitrogen storage tubes – 6000 labels and resin ribbon roll



HL0288 \$63.18
Assorted colour dots for 0.5mL to 1.5mL tubes. 5000 x dots



HT0262 \$66.39
Lab labelling tapes 16 rolls of assorted colours.



HL0298 \$20.73
Labels for cryogenic storage tubes



IP0889
\$779.00
Compact label printer with free software



Clean your lab glassware and apparatus with **QUANTUMCLEAN Concentrate**



Proudly manufactured
in Australia by
Rowe Scientific Pty Ltd



Quantumclean. For total cleaning (including radioactive decontamination) of all laboratory and industrial glassware, plastics, ceramics and stainless steel.

- **Ideal for use in ultrasonic cleaners**
- **Phosphate free**
- **Biodegradeable**

USE: Supplied as a liquid concentrate for dilution with water (approx. 2-5% is adequate for most laboratory applications). Ideal for manual cleaning or for use in an ultrasonic cleaner.



Royal Flying Doctor Service

**PROUDLY SUPPORTING THE
ROYAL FLYING DOCTOR SERVICE OF AUSTRALIA**

Your purchases of 1, 5 and 20L Quantumclean and 500mL Rowe pH buffer in the unique twin neck bottles will result in \$2 per item being sent to the RFDS to support this great Australian institution.

WIDE MOUTH BOTTLES

CODE	DESCRIPTION		PRICE
PB2319	Bottle wide mouth Boston 250mL HDPE with Polypropylene screw cap (584214)	HDPE EA	\$2.00
PB3319	Bottle wide mouth Boston 500mL Polypropylene with Polypropylene screw cap (582215)	PP EA	\$1.90
PB3390	Bottle wide mouth Boston 1L Polypropylene with Polypropylene screw cap (582216)	PP EA	\$3.90
PB0036	Bottle wide mouth Boston 2L Polypropylene with Polypropylene screw cap (t582260)	PP EA	\$11.00





A 100% Australian owned
company since 1987.

ROWE SCIENTIFIC PTY LTD

LABLIKES

ISSUE 11

Are you hitting your ATP targets?

For total confidence in your PASS and FAIL results,
look no further than Charm Sciences novaLUM IIX

July through to September, contact your local Rowe Scientific Account Manager to discuss trading up your existing luminometer

- SAVE - 15% DISCOUNT OFF NOVALUM IIX when you TRADE-UP*
- 10% off list price for 12 months forward order of PocketSwab Plus
- Real time PASS or FAIL result in 5 seconds
- Re-Test feature for FAIL results
- Photomultiplier Tube for superior light detection
- PocketSwab Surface swab – room temperature storage
- AllerGiene swab – indicator for presence of allergens
- Randomization/Shuffling of sampling plan
- novaLINK 5 software included free of charge
- Your site-specific sampling plan pre-loaded
- Lease-to-own offers available (subject to swab usage)
- All novaLUMs are programmed for Fast Alkaline Phosphatase Test (Dairy Food Safety-approved throughout Australia)



IM0216

**Limited
TRADE-UP
Offer**

GC0163

HS0169

*Conditions apply – contact us for more information and to book your no obligation free trial



FOR THE DAIRY INDUSTRY...

Do you require instrumentation and consumables for testing?

We can supply your requirements for the following applications:



- Farm gate testing for antibiotics
- Hygiene testing
- Laboratory testing
 - Cheese manufacture
 - Yoghurt manufacture
 - Liquid milk processing
 - Milk powder manufacture


Items, Consumables and Instrumentation

For testing cheese, milk and yoghurt samples. See below for an example of what we offer:



- Milk dippers
- Cheese triers
- Milk butyrometers
- Centrifuge rotors
- Cryoscope standards

- Centrifuges
- Milk sediment tester
- Jolting volumeters
- Cryoscopes





How to order online at your convenience using the Rowe Scientific Web Store and its easy to use features.

As a guest, when you arrive at the Rowe Scientific Web Store, you can search for products and populate a trolley. To obtain pricing, please request a quote on your trolley of products. Clients with a login account will see pricing directly, and will be able to place orders over the Rowe Scientific Web Store.

How to Log In

If you have a login account, enter your username and password in the login panel at the top right hand corner. This will log you into your personalised Web Store account and show your special pricing.

Your Welcome Panel

When you have successfully logged in, you will see your welcome panel to the right. Here you will find some useful features that are explained on this page.

My Backorders

My backorders is a searchable list of your current back orders with Rowe Scientific. This data is updated overnight to keep you fully informed.

Quick Order

YOUR TEN MOST PURCHASED ONLINE ITEMS
Enter quantity on any of the items in the below lists then hit the 'Add all to Trolley' button at the bottom of the page.

Items	Qty	UOM
IC19548 CONDUCTIVITY SOLUTION (0.01M KCl) (1.413mS/cm @ 25C) (*R054+)	\$20.00	1L
IC1975 CONDUCTIVITY SOLUTION (0.5M KCl) (58.670mS/cm @ 25C) (*R054+)	\$20.00	1L
FF4126 FILTER PAPER GLASS FIBRE 4.7cm WHATMAN GF/C (1822047)	\$52.41	100

Quick Order allows you place orders from a list of routinely ordered products. There are two parts to your quick order page. Your top ten recent purchases and a manual list made up of items you have added to your quick order list. These are grouped by product type e.g. Glassware for glass beakers etc.

To add a product to your quick order list, click on the add to quicklist button on the product listing. If the button does not appear on an item, it is already in your quick order list.

ADD TO QUICKLIST

To order items directly by quick order list, simply enter the quantity required for each item being ordered then select the add to trolley button. You will see these items appear in your 'trolley' panel and the total value of your order.

You can continue searching and shopping for products or view your trolley and process your order.

Order History

Here you can see a list of recent online orders placed via our Web Store.

Please note: only orders placed online appear in the list below
view last 10 orders

Order No. 211611/ Customer Ref. Carbest

Product ID	Description
FF2227	FILTER PAPER 15.0cm SIMILAR TO G4 (MN137) (4345)
FF4048	FILTER PAPER GLASS FIBRE 11.0cm SIMILAR G4 GF/B

The Search Bar

The search bar functionality on the right allows you to string key words together as you would a Google search. e.g. 'beaker 250mL' the results will list those items found which best match the exact criteria to the top of the results page. When you use the product filter buttons a second inner search bar appears in the search results page allowing you to drill down from the results using key words.

Search all stock...

Search by YOUR part number

We can upload your part numbers (where supplied by you) for contract lines allowing you to also search by your part number for easy and accurate ordering of what you need. Your part number will appear in red in brackets directly after the Rowe code.

FF2900 (10617517) FI
(MN1672) (472032)
Unit:100 Local Stock:64
Est. Lead days if no local or

See your agreed pricing

You will see your agreed contract and special pricing. This is indicated by the words **Your contract/special price** above the price.

Your Special price
\$10.00
VOLUME PRICING
QTY

Availability

You can check the availability and estimated lead times on products. Faster lead times may be available on request. The national stock is the total of our stock across the country.

TRACEABLE, IN CONTAMINATION PRO
(R608+)
Unit:500mL Local Stock:463 National Stock:76
Est. Lead days if no local or national stock:5
SSS C.O.F.A. ADD TO QUICKLIST

Account manager use.

www.rowe.com.au/rowe/KoehlerFlashPoint

USERNAME:

PASSWORD:

If you don't have a Rowe Scientific login, go to www.rowe.com.au and click the register tab to request one.



Discover more at:

www.rowe.com.au

Online, by phone or face to face, we give you the choice.



New South Wales & ACT

Ph: (02) 9603 1205
rowensw@rowe.com.au

Queensland

Ph: (07) 3376 9411
roweqld@rowe.com.au

South Australia & NT

Ph: (08) 8186 0523
rowesa@rowe.com.au

Victoria & Tasmania

Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia

Ph: (08) 9302 1911
rowewa@rowe.com.au



*A 100% Australian
owned company
since 1987.*



REF:547 Prices do not include GST. We reserve the right to change specifications, details and descriptions without notice. Pictures for illustrative purposes only. Discounts do not apply to service, freight and or repairs. Discount is off list price, does not apply to existing discounted pricing. Special pricing only while stocks last. X:\MARKETING\ Lablikes.Issue9

