



ROWE SCIENTIFIC PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



BagMixer® range By Interscience now available at Rowe Scientific.



BagMixer® 400

IS4200



**ONLY
\$3200**

Save \$270

The BagMixer® lab blender ensures optimal extraction by the action of its two paddles. Adapted to any sample (food, agricultural, chemical, medical, environmental) to be analysed between 50 mL and 400 mL, it ensures homogenisation in 30 to 60 seconds.

QUICK HOMOGENISING OF SAMPLES

- Optimal bacterial extraction
- Only 60 seconds to homogenise a sample

USER-FRIENDLY

- Compact
- Adjustable blending time (30 - 210 secs. or open running)
- Q-Tight® closing: no more leaks

EASY CLEANING

- Fully opening door (270°)
- All stainless steel

BagMixer® 400CC

IS4202



**ONLY
\$4990**

Save \$1500

GREAT DEAL!

BagMixer® 400 CC® is an easy-to-use and powerful lab blender. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses avoiding all risks of cross-contamination.

QUICK HOMOGENISING OF SAMPLES

- Optimal bacterial extraction
- Only 60 seconds to homogenise a sample

USER-FRIENDLY

- Compact
- Multifunctional screen: easy control
- Side-by-side paddle stop: easy insertion of the bag
- Q-Tight® closing: no more leaks

EASY CLEANING

- Exclusive Click & Clean® removable paddles (Patented)
- Fully opening window door (270°)
- Security drip tray
- All stainless steel

ACCURACY

- Adjustable blending speed (4 - 6 - 8 - 10 strokes per second)
- Adjustable blending time (1 sec - 1 hour or open-running)
- Adjustable blending power : adapted to any kind of sample (Patented)

www.rowe.com.au

**Special price for a limited time only.
While current stocks last.**



South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au