



ROWE SCIENTIFIC
PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



Transferpette® S

CENTRALLY LOCATED!
Large pipetting button allows uniform, smooth movement of the piston.

SECURE!
A brake helps protect against accidental volume changes.

CLEARLY ACCURATE!
4-digit volume display for precise volumes. Always remains visible during pipetting.

COLOUR-CODED
Colour indicates volume range.

SLIM
If the tip fits in the vessel, so does the shaft - no retrofitting necessary.

Maximum compatibility!

The tip cone design offers excellent fit with pipette tips from BRAND and broadset compatibility with tips from MOST OTHER manufacturers.



STARTER KIT DEALS

Brand Transferpette S Starter Kits	DESCRIPTION	ROWE CODE	LIST	NOW
Micro Kit. 3 x variable pipettes	Pipettor set Transferpette S, micro starter kit, consists of: 3 pipettes (2,5uL, 10uL & 100uL variable type), 3 x tipbox filled, 3 shelf/rack mounts	DP0073	820.00	738.00
Midi Kit. 3 x variable pipettes	Pipettor set Transferpette S, midi starter kit, consists of: 3 pipettes (20uL, 200uL & 1000uL variable type), 3 x tipbox filled, 3 shelf/rack mounts	DP0074	776.00	698.00
Macro Kit. 3 x variable pipettes	Pipettor set Transferpette S, macro starter kit, consists of: 3 pipettes (1000uL, 5000uL & 10000uL variable type), 3 x tipbox filled, 3 shelf/rack mounts	DP4922	820.00	738.00
Standard Kit. 3 x variable pipettes	Pipettor set Transferpette S, standard starter kit, consists of: 3 pipettes (10uL, 100uL & 1000uL variable type), 3 x tipbox filled, 3 shelf/rack mounts	DP0075	794.40	714.96
Pipette Package. 5 x variable pipettes	Pipettor set Transferpette S, pipette package 1, consists of: 5 pipettes (2,5uL, 10uL, 100uL, 200uL & 1000uL), 3 x tipbox filled, bench-top rack for 6 Transferpette S	DP0181	1,383.80	1,245.42

Please note limited kits available.

The Transferpette® S micropipettes are available as manual single-channel and multi-channel pipettes. They are of high quality and have many features that make pipetting easy and efficient, e.g. genuine one-hand operation, short stroke, low operating forces, volume change protection and a clearly visible volume display. They are service friendly and can be adjusted without special tools with the Easy Calibration Technology.

The single-channel pipettes Transferpette® S cover the volume range from 0.1 µl to 10 ml. The multi-channel pipettes Transferpette® 8 and -12 are available with volumes from 0.5 µl up to 300 µl.

All Transferpette® S models are completely autoclavable at 121°C, DE-M-compliant according to German Measurement and Calibration Regulation and CE-IVD-compliant according to the IVD Directive 98/79 EC. Suitable for temperatures from +15 °C to +40 °C. Vapor pressure up to 500 mbar. Viscosity: 260 mPas (260 cps).

order today - **www.rowe.com.au**

South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au